

# virotype<sup>®</sup> BTV pan/8 RT-PCR Kit Handbook

For detection of RNA from *Bluetongue Virus*  
and European serotype BTV-8

---

Licensed in accordance with § 11 (2) of the German Animal Health Act  
MA No.: FLI-B 539



24 reactions (cat. no. VT280443)



96 reactions (cat. no. VT280445)



INDICAL BIOSCIENCE GmbH, Deutscher Platz 5b,  
04103 Leipzig, Germany

# Contents

Kit contents .....	3
Intended use .....	3
Symbols .....	4
Quality control .....	4
Storage .....	5
Safety information .....	5
Introduction .....	6
Principle .....	6
RNA extraction .....	7
Equipment and reagents to be supplied by user .....	9
Important notes .....	10
General precautions .....	10
Protocol: Real-time RT-PCR for detection of RNA from <i>Bluetongue Virus</i> and BTV-8 .....	12
Important points before starting .....	12
Things to do before starting .....	12
Procedure when using 96-well plate real-time cycler .....	12
Procedure when using Rotor-Gene® Q or similar thermal cyclers ...	15
Data analysis and interpretation .....	18
Interpretation of results .....	18
Change index .....	23

## Kit contents

<b>virotype BTV pan/8 RT-PCR Kit</b>	<b>(24)</b>	<b>(96)</b>
<b>Cat. no.</b>	<b>VT280443</b>	<b>VT280445</b>
<b>Number of reactions</b>	<b>24</b>	<b>96</b>
Master Mix (tube with orange cap) includes enzymes, primers, and probes	1 x 500 µl	2 x 980 µl
Positive Control (tube with red cap)	1 x 25 µl	1 x 150 µl
Negative Control (tube with blue cap)	1 x 25 µl	1 x 150 µl
Handbook	1	1

## Intended use

The virotype BTV pan/8 RT-PCR Kit is intended for the detection of *Bluetongue Virus* RNA in ruminant whole blood (preferred with anticoagulants, for example EDTA-blood) and tissue samples (spleen, lymph nodes) from cattle, sheep, and goats.

The kit is approved by the Friedrich-Loeffler-Institut and licensed in accordance with § 11 (2) of the German Animal Health Act (FLI-B 539) for use in Germany for veterinary diagnostic procedures.

**For veterinary use only.**

# Symbols



Legal manufacturer



Lot number



Use by date



Temperature limitations for storage



Handbook



Catalog number



Material number



Protect from light



For cattle samples

# Quality control

In accordance with INDICAL's ISO-certified Quality Management System, each lot of virotype BTV pan/8 RT-PCR Kit is tested against predetermined specifications to ensure consistent product quality.

## Storage

The components of the virotype BTV pan/8 RT-PCR Kit should be stored at -30°C to -15°C and are stable until the expiration date stated on the label. Avoid repeated thawing and freezing (>2x), as this may reduce assay sensitivity. Freeze the components in aliquots if they will only be used intermittently.

## Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available from your local sales representative or by Email request under **compliance@indical.com**.

All sample residues and objects that have come into contact with samples must be decontaminated or disposed of as potentially infectious material.

# Introduction

Bluetongue is an infectious, non-contagious disease of ruminants. The agent is the *Bluetongue Virus* (BTV), a double-stranded RNA virus of the genus *Orbivirus* of the family *Reoviridae* which includes at least 25 known serotypes. BTV is widely distributed around the world. Sheep, cattle, and goats are mainly affected by the disease. Clear clinical signs are usually seen only in sheep. In severe cases the tongue may show intense hyperemia and become cyanotic (Bluetongue).

BTV serotype 8 is of epidemiological importance in Central Europe and cause of recent major Bluetongue Disease outbreaks. The virus is transmitted by certain midges of the genus *Culicoides*. Furthermore, the virus can be spread by contaminated needles and surgery equipment.

## Principle

Polymerase chain reaction (PCR) is based on the amplification of specific regions of the pathogen genome. In real-time RT-PCR, the amplified product is detected using fluorescent dyes. These are usually linked to oligonucleotide probes that bind specifically to the amplified product. Monitoring the fluorescence intensities during the PCR run (e.g., in real time) allows the detection of the accumulating product without the need to re-open the reaction tubes afterward.

The virotype BTV pan/8 RT-PCR Kit contains all of the necessary reagents for the detection of BTV RNA, including a Positive and Negative Control. With this kit, both reverse transcription and PCR are performed in one reaction tube, reducing the risk of contamination.

The virotype BTV pan/8 RT-PCR Kit uses three specific primer/probe combinations:

- FAM™ fluorescence for RNA of the 25 known BTV serotypes
- Cy®5 fluorescence for RNA of European serotype 8 (BTV-8)
- HEX™ fluorescence for the internal control ( $\beta$ -actin mRNA present within the sample)

The Positive Control contains BTV-8 RNA and allows the control of the denaturation step since the successful denaturation of the viral double-stranded RNA is a prerequisite for amplification.

## RNA extraction

virotype BTV pan/8 can be used for the detection of BTV RNA from ruminant whole blood (preferred with anticoagulants, e.g. EDTA-blood) and tissue samples (spleen, lymph nodes).

Due to the high sensitivity of the test, pools of up to 10 individual blood samples may be analyzed. However, the optimum pool size depends on the regional prevalence for BTV.

Prior to real-time RT-PCR, viral RNA must be extracted from the starting material. INDICAL offers a range of validated kits for the extraction of RNA from animal samples.

### Extraction based on magnetic beads:

- MagAttract® 96 cador® Pathogen Kit

Extraction based on spin columns:

- QIAamp® cador® Pathogen Mini Kit
- cador Pathogen 96 QIAcube® HT Kit

If real-time RT-PCR is not performed immediately after extraction, store the RNA at -20°C or at -70°C for longer storage.

For further information on automated and manual extraction of BTV RNA from different sample types, refer to the respective handbook or contact INDICAL Support and **support@indical.com**.

## Equipment and reagents to be supplied by user

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, consult the appropriate safety data sheets (SDSs), available from the product supplier.

- Pipets
- Nuclease-free, aerosol-resistant pipet tips with filters
- Sterile 1.5 ml Eppendorf® tubes
- Nuclease-free (RNase/DNase-free) consumables. Special care should be taken to avoid nuclease contamination of all reagents and consumables used to set up PCR for sensitive identification of viral nucleic acids
- Cooling device or ice
- Benchtop centrifuge with rotor for 1.5 ml tubes
- Real-time cycler with appropriate fluorescent channels
- Appropriate software for chosen real-time cycler
- Appropriate strip tubes and caps or 96-well optical microplate with optical sealing film or cover for chosen real-time cycler

# Important notes

## General precautions

The user should always pay attention to the following:

- Use nuclease-free pipet tips with filters.
- Store and extract positive materials (specimens, positive controls and amplicons) separately from all other reagents, and add them to the reaction mix in a spatially separated facility.
- Thaw all components on ice before starting as assay.
- When thawed, mix the components by inverting and centrifuge briefly.
- Do not use components of the test kit past the expiration date.
- Keep samples and controls on ice or in a cooling block during the setup of reactions.

## Negative control

At least one negative control reaction should be included in each PCR run. This enables assessment of contamination in the reaction.

## Positive control

When performing PCR on unknown samples, it is recommended to perform a positive control reaction in the PCR run, containing a sample that is known to include the targeted viral RNA. A positive control serves to prove the functionality of the pathogen assay, for example, the correct setup of the reaction mix. Use 5 µl of the Positive Control provided with the virotype BTV pan/8 RT-PCR Kit to test for successful amplification of the target.

## Extraction and amplification control

For increased process safety and convenience, an extraction and amplification control assay is included in the form of an additional primer/probe set that detects a housekeeping gene present within the sample. This allows both extraction and amplification to be monitored.

# Protocol: Real-time RT-PCR for detection of RNA from *Bluetongue Virus* and BTV-8

## Important points before starting

- Please read „Important notes“ on page 10 before starting.
- Include at least one positive control (Positive Control) and one negative control (Negative Control) per PCR run.
- Before beginning the procedure, read through the protocol and ensure that you are familiar with the operation of the chosen real-time PCR cycler.
- RNA is unstable. Perform the protocol without interruption.

## Things to do before starting

- Thaw all reagents on ice and protect from light.
- Maintain reagents on ice during PCR setup.
- Before use, spin the reagents briefly.

## Procedure when using 96-well plate real-time cycler

1. Pipet 5 µl of RNA samples, Positive Control, and Negative Control into individual reaction tubes. Cover the reaction tubes (e.g., with PCR sealing foil).

Include positive and negative control reactions.

Positive control: Use 5 µl of the positive control (Positive Control) instead of sample RNA.

- Negative control: Use 5 µl of the negative control (Negative Control) instead of sample RNA
- Denature the samples for 5 min at 98°C in a 96-well plate standard cycler with a heated lid.
  - Immediately cool down on ice water or liquid nitrogen for at least 20 s. Then store on ice or in a cooling device.
  - Pipet 20 µl of the Master Mix into each reaction tube. Thus the final volume of a test is 25 µl (Table 1).

Table 1. Preparation of reaction mix

<b>Component</b>	<b>Volume</b>
Master Mix	20 µl
Sample	5 µl
<b>Total volume</b>	<b>25 µl</b>

- Close the reaction tubes with the corresponding caps.
- Set the filters for the reporter dyes in the software of your thermal cycler according to Table 2.

Table 2. Filter settings for the reporter

<b>Pathogen/ Internal Control</b>	<b>Reporter</b>
BTV	FAM
BTV-8	Cy5
Internal Control	HEX/ JOE™ <sup>1</sup>
Passive reference <sup>2</sup>	ROX™

<sup>1</sup> Use the option appropriate for your thermal cycler.

<sup>2</sup> Internal reference for use on ABI PRISM® Sequence Detection Systems by Applied Biosystems®

7. Run the real-time RT-PCR protocol according to Table 3 if running only the virotype BTV pan/8 RT-PCR Kit.

Table 3. Real-time RT-PCR protocol for BTV pan/8

<b>Step</b>	<b>Temperature</b>	<b>Time</b>	<b>Number of cycles</b>
<b>Reverse Transcription</b>	50°C	10 min	1
<b>Initial Activation</b>	95°C	10 min	1
<b>2-step cycling</b>			
Denaturation	95°C	15 s	40
Annealing/Extension*	60°C	60 s	

\* Fluorescence data collection.

8. Run the real-time RT-PCR protocol according to Table 4 if running other virotype assays simultaneously (e.g., virotype BTV pan/4, virotype BVDV and/ or virotype SBV).

Table 4. Real-time RT-PCR protocol for simultaneous assays

<b>Step</b>	<b>Temperature</b>	<b>Time</b>	<b>Number of cycles</b>
<b>Reverse Transcription</b>	50°C	20 min	1
<b>Initial Activation</b>	95°C	15 min	1
<b>3-step cycling</b>			
Denaturation	95°C	30 s	40
Annealing*	57°C	45 s	
Extension	68°C	45 s	

## Procedure when using Rotor-Gene® Q or similar thermal cyclers

1. Pipet at least 7 µl of RNA samples or Positive Control into individual 0.2 ml PCR reaction tubes. Cover the reaction tubes (e.g., with PCR sealing foil).

Include positive control reactions.

Positive control: Use at least 7 µl of the positive control (Positive Control) instead of sample RNA.

2. Denature the samples for 5 min at 98°C in a 96-well plate standard cyler with a heated lid.
3. Immediately cool down on ice water or liquid nitrogen for at least 20 s. Then store on ice or in a cooling device.
4. Pipet 5 µl of RNA samples, Positive Control, and Negative Control into individual strip tubes (0.1 ml), which are suitable for the Rotor-Gene Q or chosen thermal cyler.
5. Add 20 µl of the Master Mix into each reaction tube. Thus the final volume of a test is 25 µl (Table 5).

Table 5. Preparation of reaction mix

<b>Component</b>	<b>Volume</b>
Master Mix	20 µl
Sample	5 µl
<b>Total volume</b>	<b>25 µl</b>

6. Close the reaction tubes with the corresponding caps.
7. Set the filters for the reporter dyes in the software of your thermal cyler according to Table 6.

Table 6. Filter settings for the reporter using the Rotor-Gene Q

<b>Pathogen/ Internal Control</b>	<b>Reporter</b>
BTV	green/ FAM
BTV-8	red/ Cy5
Internal Control	yellow/ HEX

- Run the real-time RT-PCR protocol according to Table 7 if running only the virotype BTV pan/8 RT-PCR Kit.

Table 7. Real-time RT-PCR protocol for BTV pan/8

<b>Step</b>	<b>Temperature</b>	<b>Time</b>	<b>Number of cycles</b>
<b>Reverse Transcription</b>	50°C	10 min	1
<b>Initial Activation</b>	95°C	10 min	1
<b>2-step cycling</b>			
Denaturation	95°C	15 s	40
Annealing/Extension*	60°C	60 s	

\* Fluorescence data collection.

- Run the real-time RT-PCR protocol according to Table 8 if running other virotype assays simultaneously (e.g., virotype BTV pan/4, virotype BVDV and/ or virotype SBV).

Table 8. Real-time RT-PCR protocol for simultaneous assays

<b>Step</b>	<b>Temperature</b>	<b>Time</b>	<b>Number of cycles</b>
<b>Reverse Transcription</b>	50°C	20 min	1
<b>Initial Activation</b>	95°C	15 min	1
<b>3-step cycling</b>			
Denaturation	95°C	30 s	
Annealing*	57°C	45 s	40
Extension	68°C	45 s	

\* Fluorescence data collection.

# Data analysis and interpretation

## Interpretation of results

For the assay to be valid the Positive Control must give a signal in the FAM, Cy5 and HEX channels with a  $C_T^1 < 35$ . If no FAM and no Cy5 signals of the Positive Controls are measured the denaturation and cooling steps were insufficient and the testing should be repeated. The Negative Control must give no signal.

The following results are possible if working with unknown samples. The possible sample results are also summarized in Table 9 on page 20.

**The sample is positive for BTV and BTV-8, and the assay is valid, if the following criteria are met:**

- The sample yields a signal in the FAM, Cy5 and HEX channel.
- The Positive Control yields a signal in all channels.
- The Negative Control yields no signal in any of the channels.

Note that very high concentrations of BTV RNA in the sample may lead to a reduced HEX signal or no HEX signal due to competition with the internal control.

**The sample is positive for BTV and negative for BTV-8, and the assay is valid, if the following criteria are met:**

- The sample yields a signal in the FAM and HEX channel, but not in the Cy5 channel.
- The Positive Control yields a signal in all channels.

---

<sup>1</sup> Threshold cycle ( $C_T$ ) — cycle at which the amplification plot crosses the threshold, i.e., there is the first clearly detectable increase in fluorescence

- The Negative Control yields no signal in any of the channels.

Note that very high concentrations of BTV RNA in the sample may lead to a reduced HEX signal or no HEX signal due to competition with the internal control.

**The sample is negative for both BTV and BTV-8, and the assay is valid, if the following criteria are met:**

- The sample yields a signal only in the HEX channel.
- The Positive Control yields a signal in all channels.
- The Negative Control yields no signal in any of the channels.

A positive HEX signal rules out the possibility of PCR inhibition and/ or incorrect RNA extraction as the internal control is amplified.

**The sample results are inconclusive, and the assay is invalid, if the following occurs:**

- The sample yields no signal in any of the fluorescence channels.

The PCR was inhibited or the sample extraction was incorrect. It is recommended to retest the respective individual samples in nuclease free water (e.g., diluted 1:5), to repeat the RNA extraction, or repeat the whole test procedure starting with new sample material.

Check that there is a fluorescence signal in the all channels for the positive control reaction (Positive Control). Absence of a signal for the Positive Control indicates an error, which could be due to incorrect RNA denaturation, RNA extraction failure,, or incorrect cycling conditions.

Repeat RNA extraction or repeat the whole procedure starting with new sample material.

Table 9. Results interpretation table\*

<b>FAM</b>	<b>Cy5</b>	<b>HEX</b>	<b>Sample result</b>
			Positive for:
X		(X)	BTV
X	X	(X)	BTV <u>and</u> BTV-8
		X	negative
			inconclusive

\* Interpretation of sample results can be determined provided positive and negative control reactions are performed. The Positive Control must yield a signal in the FAM, Cy5 and HEX channels. The Negative Control must yield no signal. For a complete explanation of possible sample results please refer to “Data analysis and interpretation” on page 18.

INDICAL offers a range of ELISA kits and real-time PCR and real-time RT-PCR kits for the detection of animal pathogens.

Visit **[www.indical.com](http://www.indical.com)** for more information about bactotype, cadon, cattletype, flocktype, pigtype and virotype products.

For up-to-date licensing information and product-specific disclaimers, see the respective INDICAL kit handbook or user manual.

## Notes

## Limited License Agreement for virotype BTV pan/8 RT-PCR Kit

Use of this product signifies the agreement of any purchaser or user of the product to the following terms:

1. The product may be used solely in accordance with the protocols provided with the product and this handbook and for use with components contained in the kit only. INDICAL grants no license under any of its intellectual property to use or incorporate the enclosed components of this kit with any components not included within this kit except as described in the protocols provided with the product, this handbook, and additional protocols available at [www.indical.com](http://www.indical.com). Some of these additional protocols have been provided by INDICAL users for INDICAL users. These protocols have not been thoroughly tested or optimized by INDICAL. INDICAL neither guarantees them nor warrants that they do not infringe the rights of third-parties.
2. Other than expressly stated licenses, INDICAL makes no warranty that this kit and/or its use(s) do not infringe the rights of third-parties.
3. This kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. INDICAL specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. INDICAL may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the kit and/or its components.

For updated license terms, see [www.indical.com](http://www.indical.com).

**Trademarks:** bactotype<sup>®</sup>, cador<sup>®</sup>, cattletype<sup>®</sup>, flocktype<sup>®</sup>, pigtype<sup>®</sup>, virotype<sup>®</sup> (INDICAL BIOSCIENCE GmbH); MagAttract<sup>®</sup>, QIAamp<sup>®</sup>, QIAcube<sup>®</sup>, Rotor-Gene<sup>®</sup> (QIAGEN GmbH); ABI PRISM<sup>®</sup> 7500 (Applied Biosystems); FAM<sup>™</sup>, HEX<sup>™</sup>, JOE<sup>™</sup>, ROX<sup>™</sup> (Life Technologies Corporation); Cy<sup>®</sup>5 (GE Healthcare); Eppendorf<sup>®</sup> (Eppendorf-Netheler-Hinz GmbH). Licensed probes manufactured by Integrated DNA Technologies, Inc. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

HB-1598-EN-005 © 2018 INDICAL BIOSCIENCE GmbH, all rights reserved.

## Change index

Handbook	Version	Change
HB-1598-EN-005	August 2018	INDICAL design

**INDICAL**  
BIOSCIENCE

Ordering: [www.indical.com/contact](http://www.indical.com/contact)  
Technical Support: [support@indical.com](mailto:support@indical.com)  
Website: [www.indical.com](http://www.indical.com)