

virotype[®] BVDV Positive Control Product Sheet



60 reactions (cat. no. PC280377)



INDICAL BIOSCIENCE GmbH, Deutscher Platz 5b,
04103 Leipzig, Germany

Contents

virotype BVDV Positive Control	(60)
Cat. no.	PC280377
Number of reactions	60
BVDV Positive Control	2 x 150 µl

Intended use

The virotype BVDV Positive Control contains synthetically derived RNA of *Bovine Viral Diarrhea Virus*. This product can be used as a positive control in real-time RT-PCR.

For research use only. This product is not intended for the diagnosis, prevention or treatment of a disease.

Preparation of PCR setup

Before use, thaw the reagent and spin down briefly.

The Positive Control can be used with a variety of PCR chemistries.

Product Name	Cat No	Volume recommendation *
virotype Mix +IC(JOE)-RNA	MX299935	8 µl
IndiMix JOE	MX299945	8 µl
virotype RT/Mix +IC-RNA	MX280377	5 µl

* per 25 µl PCR reaction

Quality control

Each lot of this product was manufactured and released in accordance with INDICAL's ISO-certified Quality Management System.

Storage and handling

Upon receipt, store the material at -30°C to -15°C, away from any sources of contaminating DNA or RNA, especially amplified DNA products and material used for PCR master mixes. Use aerosol barrier pipet tips for pipetting. Do not use after the expiration date printed on the label.

Avoid repeated thawing and freezing because this may reduce performance. If used intermittently, aliquot the material and freeze.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available from your regional sales manager or at **compliance@indical.com**.

All sample residues and objects that have come into contact with samples must be decontaminated or disposed of as potentially infectious material.

Warning and precautions

Dispose of all materials in accordance with national and local regulations.

Trademarks: virotype® (INDICAL BIOSCIENCE GmbH). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. For updated license terms, see www.indical.com.

HB-2501-EN-002 © 2018-2022 INDICAL BIOSCIENCE GmbH, all rights reserved.

Change index

Product Sheet	Version	Change
HB-2501-EN-002	June 2022	Revision Chapter PCR Setup
HB-2501-EN-001	October 2018	Product launch