

IndiSpin Pathogen Kit

Quick Start Protocol

For manual use.

The IndiSpin Pathogen Kits (cat. no. SP54104 and SP54106) can be stored at room temperature (15-25°C). For expiry date information, read the label on the kit box.

Further information and support

- IndiSpin Pathogen Kit Handbook: www.indical.com/handbooks
- Technical assistance: support@indical.com

Important notes before starting

- Dissolve carrier RNA in Buffer AVE as indicated on the tube.
- Add isopropanol (100%) to Buffer ACB and ethanol (96-100%) to Buffers AW1 and AW2 before use. See the respective bottle labels for volumes.
- Carry out all centrifugation steps at room temperature in a conventional table-top microcentrifuge.
- Equilibrate buffers to room temperature (15-25°C).

For up-to-date licensing information and product-specific disclaimers, see the respective INDICAL kit handbook or user manual. Trademarks: cadon[®], IndiSpin[®] and INDICAL[®] (INDICAL BIOSCIENCE).

Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are not to be considered unprotected by law. SAP: 1116318; HB-2517-EN-001 10/2018 © 2018 INDICAL BIOSCIENCE GmbH, all rights reserved.

Procedure

1. Pipet 20 μ l Proteinase K into a 2 ml microcentrifuge tube (not provided).
2. Add 200 μ l fluid sample to the Proteinase K.
Note: If your sample volume is less than 200 μ l, bring it to 200 μ l by adding PBS.
3. Add 100 μ l Buffer VXL. Close the cap and mix by pulse vortexing.
Note: For cell-free samples, ensure that 1 μ g Carrier RNA is added per 100 μ l Buffer VXL.
4. Incubate at 20-25°C for 15 min.
5. Briefly centrifuge the tube to remove drops from the inside of the lid.
6. Add 350 μ l Buffer ACB to the sample and mix thoroughly by pulse-vortexing.
7. Briefly centrifuge the tube to remove drops from the inside of the lid.
8. Transfer the lysate to a spin column placed in a 2 ml collection tube. Close the cap and centrifuge at 6000 x g (8000 rpm) for 1 min. Transfer the spin column to a clean 2 ml collection tube and discard the collection tube containing the filtrate.
9. Add 600 μ l Buffer AW1 and centrifuge at 6000 x g (8000 rpm) for 1 min. Transfer the spin column to a clean 2 ml collection tube and discard the collection tube containing the filtrate.
10. Add 600 μ l Buffer AW2 and centrifuge at 6000 x g (8000 rpm) for 1 min. Transfer the spin column to a clean 2 ml collection tube and discard the collection tube containing the filtrate.
11. Centrifuge the spin column at 20,000 x g (14,000 rpm) for 2 min.
12. Place the spin column in a clean 1.5 ml microcentrifuge tube (not provided) and discard the collection tube containing the filtrate. Add 50-150 μ l Buffer AVE to the center of the membrane, close the cap and incubate at room temperature for 1 min.
13. Centrifuge at 20,000 x g (14,000 rpm) for 1 min.