SVANOVIR® BLV gp51-Ab

Bovine Leukaemia Virus gp51

Antibody Test

Screening format

Bovines Leukämie-Virus gp51

Antikörper Test

Monophasisch

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Contents	Art. No. SV-104886
Microtitre plate Microtitre plates (96 wells) coated with non-infectious BLV gp51 antigen (sealed and stored dry)	10 (Strips) 12 x 8
Conjugate Ready-to-use (horseradish peroxidase conjugated anti-bovine IgG monoclonal antibodies)	1 x 120 mL
PBS-Tween Solution 20 x concentrate	3 x 125 mL
Sample Dilution Buffer Ready-to-use	1 x 120 mL
Substrate Solution (Tetramethylbenzidine in substrate buffer containing H ₂ O ₂) - STORE IN THE DARK!	1 x 100 mL
Stop Solution Contains sulphuric acid (2M) - DANGER!	 2 x 25 mL
A. Positive Control Serum - Contains preservatives	1 x 0.5 mL
B. Negative Control Serum - Contains preservatives	1 x 0.5 mL

This manual covers the following
 SVANOVIR® BLV gp51-Ab kit:
 Article number SV-104886

Bovine Leukaemia Virus gp51

Antibody Test

Screening format

Name and Application

The SVANOVIR® BLV gp51-Ab is an Enzyme Linked Immunosorbent Assay (ELISA) for the detection of antibodies against BLV gp51 in bovine serum, plasma and milk samples; individualsamples, pools and bulk tank milk.

General information

Bovine leukaemia virus (BLV) is an exogenous oncogenic retrovirus infecting B lymphocytes. The infection may lead to persistent lymphocytosis and, after lifelong latency, *i.e.* in some adult cattle, to development of tumours (lymphosarcomas) with associated symptoms^{1,2}. The spread does often start from infected cattle introduced into herds and then take enzootic proportions. This has given the disease the name 'enzootic bovine leucosis' - in contrast to sporadic or juvenile bovine leucosis of unknown etiology. Since BLV is bound to lymphocytes, the exact mode of transmission by direct contact is unknown. In some countries transmission by horseflies is anticipated³. Transplacental transmission and transmission via milk are considered uncommon. On the contrary, iatrogenic transmission occurs easily since as little as 0.1 µL blood from an infected cow did transmit infection⁴. There is an antibody response to infection, and antibodies are found in both serum and milk. The routine diagnosis of BLV infection is thus based on antigen-antibody reaction⁵. Eradication by culling, if planned by a country with a high BLV-infection rate, is a formidable task. The control program usually starts from the results of a country-wide screening for antibodies to BLV, particularly in dairy herds.

Principle

The kit procedure is based on a solid phase indirect ELISA. A monoclonal antibody against gp51 is adsorbed to microtiter plates or strips. This antibody captures gp51 antigen. BLV antibodies (if present in the test sample) bind to the antigen in the well. HRP conjugate added subsequently forms a complex with the BLV antibodies. Unbound material is removed by rinsing before the addition of a substrate solution. Subsequently a blue colour develops which is due to the conversion of the substrate by the conjugate. A positive result is indicated by development of a blue colour. The reaction is stopped by addition of the stop solution; the colour changes to yellow. The result can be read visually or by a microplate photometer, where the optical density (OD) is measured at 450 nm.

Materials needed but not provided

1. Precision pipettes
2. Disposable pipette tips
3. Distilled, deionised or any similar high quality water
4. Wash bottle, multichannel pipetors or plate washer
5. Container: 1 to 2 litres for PBS-Tween
6. Microplate photometer, 450 nm filter

Specimen information

Individual serum or pools up to 10 samples:

4 μ L of blood serum or plasma is needed for each sample well. Fresh, refrigerated or previously frozen serum or plasma may be tested.

Individual or pooled milk samples/bulk tank milk up to at least 50 animals:

100 μ L of skim milk is required for each sample well. Fresh, refrigerated, or previously frozen milk samples may be tested. It is recommended to centrifuge milk samples for 15 minutes at 2000 x g to remove the lipid layer, or leave the milk samples until the fat layer is formed on top of the sample. Pipette under the fat layer.

Preparation of reagents

PBS-Tween Buffer:

Dilute the PBS-Tween Solution 20 x concentrate 1/20 in distilled water. Prepare 500 mL per plate by adding 25 mL PBS-Tween solution to 475 mL distilled water and mix thoroughly.

N.B. Please check that there is no crystal precipitation in the bottle. If crystals are seen, please warm and shake well.

Precautions

1. Carefully read and follow all instructions.
2. Store the kit and all reagents at 2-8°C (36-46°F).
3. All reagents should equilibrate to room temperature, 18 -25°C (64-77°F) before use.
4. Handle all materials according to the Good Laboratory Practice.
5. Do not mix components or instruction manuals from different test kit batches.
6. Care should be taken to prevent contamination of kit components.
7. Do not use test kit beyond date of expiry.
8. Do not eat, drink, or smoke where specimens or kit reagents are handled.
9. Use a separate pipette tip for each sample.
10. Do not pipette by mouth.
11. Include positive and negative control serum on each plate or test strip series.
12. Use only distilled, deionised or any similar high quality water for preparation of reagents.
13. When preparing the buffers etc., measure the required volume.
14. The Stop Solution contains sulphuric acid, which is corrosive.*
15. All unused biological materials should be disposed according to the local, regional and national regulations.

Recommendations!

There is always a surplus volume for the liquid reagent. The volume mentioned on the label is the minimum obtainable.

Strips with broken seal can be stored 2-8°C (36-46°F) for up to 4 weeks.

Procedure

1. All reagents should equilibrate to room temperature 18-25°C (64-77°F) before use. Label each strip with a number.
2. Add samples
The provided negative and positive control sera are used for both serum and milk testing.

Serum Samples

- A. Add 100 µL of Sample Dilution Buffer to each well that will be used for serum samples and serum controls.
 - B. Add 4 µL of Positive Control Serum (Reagent A) and 4 µL of Negative Control Serum (Reagent B) respectively to selected wells coated with BLV gp51 viral antigen.
For confirmation purposes it is recommended to run the samples in duplicates.
 - C. Add 4 µL of serum sample to a selected well coated with BLV gp51 viral antigen. The samples can be tested in singlicates or in duplicates. However for confirmation purposes it is recommended to run the samples in duplicates.
3. Shake the plate thoroughly. Seal the plate/strip and incubate at 37°C (98.6°F) for 1 hour.
 4. Rinse the plates/strips 3 times with PBS-Tween Buffer: At each rinse cycle fill up the wells, empty the plate and tap hard to remove all remains of fluid.
 5. Add 100 µL of HRP Conjugate to each well and incubate at 37°C (98.6°F) for 1 hour.
 6. Repeat step #4.
 7. Add 100 µL Substrate Solution to each well. Incubate for 10 minutes at room temperature, 18-25°C (64-77°F). Begin timing when the first well is filled.
 8. Stop the reaction by adding 50 µL of Stop Solution to each well and mix thoroughly. Add the Stop Solution in the same order as the Substrate Solution in step #7.
 9. Measure the optical density (OD) of the controls and samples at 450 nm in a microplate photometer (use air as blank). Measure the OD within 15 minutes after the addition of Stop Solution to prevent fluctuation in OD values.

Continue at step #3

Milk Samples

- A. For addition of controls, see "serum Samples" (point A and B).
- B. Add 100 µL of milk sample to a selected well coated with BLV gp51 viral antigen. The samples can be tested in singlicates or in duplicates. However for confirmation purposes it is recommended to run the samples in duplicates.

Continue at step #3

Calculations

Calculations of results are done as described below.

Calculation of Percent Positivity Values (PP)

All OD Values for the test samples as well as the Negative Control are related to the OD value of the positive control as follows:

$$PP = \frac{OD_{\text{Sample or Negative control}}}{OD_{\text{Positive control}}} \times 100$$

Interpretation of the results

Criteria for test validity

To ensure validity, the duplicate OD values of the positive control should not differ more than 25% from the mean value of the two duplicates.

Additionally, the control values should fall within the following limits:

OD Positive Control > 1.0

PP Negative Control < 15

Should any of these criteria not be fulfilled, the test is invalid. For invalid tests, technique may be suspect and the assay should be repeated.

To confirm test sample results, each separate duplicate PP values should be equally interpreted, *i.e.* positive and negative. In case of discrepancy it is recommended to re-test the sample.

Interpretation of test sample result

Sample	PP	Interpretation
Serum	< 20 ≥ 20	Negative Positive
Milk	< 15 ≥ 15	Negative Positive

References

1. Burny, A., Bruck, C. *et al.* (1980) Bovine leukemia virus: molecular biology and epidemiology. In *Viral Oncology*. Edited by G. Klein. New York, Raven Press. 231-289.
2. Ferrer, J.,F. (1980) Bovine lymphosarcoma. *Adv. Vet. Sci. Comp. Med.* 24, 1-68.
3. Manet, G. *et al.* (1989). Natural mode of horizontal transmission of bovine leukemia virus (BLV): the potential role of tabanids (*Tabanus spp.*). *Vet. Immunol. Immunopathol.* 22, 255-263.
4. Evermann, J.F. *et al.* (1986). Transmission of bovine leucosis virus by blood inoculation. *Am. J. Vet. Res.* 47, 1885-1887.
5. Hoff-Joergensen, R. (1989). An international comparison of different laboratory tests for the diagnosis of bovine leucosis: suggestions for international standardization. *Vet. Immunol. Immunopathol.* 22, 293-297.

***DANGER: Stop solution (sulphuric acid)**

May be corrosive to metals. Causes skin irritation. Causes serious eye irritation.

Keep only in original container. Wear eye protection/ face protection. Wear protective gloves.

IN CASE OF CONTACT WITH EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/ physician. If eye irritation persists: Get medical advice/ attention.

IN CASE OF CONTACT WITH SKIN: Wash with plenty of soap and water. Take off contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention. Absorb spillage to prevent material damage.

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INDICAL BIOSCIENCE GmbH