SVANOVIR® BVDV-Ab

## **Bovine Virus Diarrhoea Virus**

### **Antibody Test**

*Screening format*

## **Bovine Virus Diarrhoea Virus**

### **Antikörper Test**

*Monophasisch*

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<b>Contents</b>	<b>Art. No. SV-104883</b>
<b>Microtitre plate</b> Microtitre plates (96 wells) coated with non-infectious BVDV antigen (sealed and stored dry)	10 (Strips) 12 x 8
<b>Conjugate</b> Ready-to-use (horseradish peroxidase conjugated anti-bovine IgG monoclonal antibodies)	1 x 120 mL
<b>PBS-Tween Solution</b> 20 x concentrate	2 x 125 mL
<b>Sample Dilution Buffer</b> Ready-to-use	1 x 120 mL
<b>Substrate Solution</b> (Tetramethylbenzidine in substrate buffer containing H <sub>2</sub> O <sub>2</sub> ) - STORE IN THE DARK!	1 x 100 mL
<b>Stop Solution</b> Contains sulphuric acid (2M) - DANGER!	 2 x 25 mL
<b>A. Positive Control Serum</b> - Contains preservatives	1 x 0.7 mL
<b>B. Negative Control Serum</b> - Contains preservatives	1 x 0.7 mL

This manual covers the following  
 SVANOVIR® BVDV-Ab kit:  
 Article number SV-104883

# Bovine Virus Diarrhoea Virus

## Antibody Test

### Screening format

#### **Name and Application**

The SVANOVIR® BVDV-Ab is an Enzyme Linked Immunosorbent Assay (ELISA) for the detection of BVDV specific antibodies in bovine serum, plasma and milk samples; individual samples, pools and bulk tank milk.

#### **General information**

The Bovine Virus Diarrhoea Virus (BVDV), a pesti-virus, is associated with a range of diseases affecting the alimentary tract of cattle. Clinical signs include pyrexia, diarrhoea and reduced milk yield; there is a high morbidity but low mortality of infected animals<sup>1</sup>. Infection of pregnant cows may result in transplacental foetal infection. Foetuses may be aborted, mummified, stillborn or born with severe anomalies. However, in many instances calves, if immunotolerant, are born unaffected but with persistent viraemia and become transmission sources of the virus<sup>2</sup>. Moreover, recurrent infection of such calves may result in mucosal disease<sup>3</sup> – a condition with low morbidity but high mortality; it is characterized by extensive erosions in the oral and gastrointestinal mucosa. BVDV has also been shown to have immunosuppressive effects which predisposes animals to infection by other microorganisms<sup>4</sup>. Economic losses can be high due to these factors which is why herds should be monitored for the presence of infection.

#### **Principle**

The kit procedure is based on a solid phase indirect ELISA. In this procedure, samples are exposed to non-infectious BVDV antigen coated wells on microtitre plates or strips. BVDV antibodies (if present in the test sample) bind to the antigen in the well. The HRP conjugate added subsequently forms a complex with these BVDV antibodies. Unbound material is removed by rinsing before the addition of a substrate solution. Subsequently a blue colour develops which is due to the conversion of the substrate by the conjugate. A positive result is indicated by development of the blue colour. The reaction is stopped by addition of stop solution; the colour changes to yellow. The result can be read visually or by a microplate photometer, where the optical density (OD) is measured at 450 nm.

## Materials needed but not provided

1. Precision pipettes
2. Disposable pipette tips
3. Distilled, deionised or any similar high quality water
4. Wash bottle, multichannel pipettor or plate washer
5. Container: 1 to 2 litres for PBS-Tween
6. Microplate photometer, 450 nm filter
7. Plate shaker

## Specimen information

### Individual serum or pools up to 10 samples:

10 µL of blood serum or plasma is needed for each sample well. Fresh, refrigerated or previously frozen serum or plasma may be tested.

### Individual or pooled milk samples/bulk tank milk up to 50 animals:

100 µL of skim milk is required for each sample well. Fresh refrigerated or previously frozen milk samples may be tested. It is recommended to centrifuge milk samples for 15 minutes at 2000 x g to remove the lipid layer, or leave the milk samples until the fat layer is formed on top of the sample. Pipette under the fat layer.

## Preparation of reagents

### PBS-Tween Buffer:

Dilute the PBS-Tween Solution 20 x concentrate 1/20 in distilled water. Prepare 500 mL per plate by adding 25 mL PBS-Tween solution to 475 mL distilled water and mix thoroughly.

**N.B.** Please check that there is no crystal precipitation in the bottle. If crystals are seen, please warm and shake well.

## Precautions

1. Carefully read and follow all instructions.
2. Store the kit and all reagents at 2-8°C (36-46°F).
3. All reagents should equilibrate to room temperature 18-25°C (64-77°F) before use.
4. Handle all materials according to the Good Laboratory Practice.
5. Do not mix components or instruction manuals from different test kit batches.
6. Care should be taken to prevent contamination of kit components.
7. Do not use test kit beyond date of expiry.
8. Do not eat, drink, or smoke where specimens or kit reagents are handled.
9. Use a separate pipette tip for each sample.
10. Do not pipette by mouth.
11. Include positive and negative control serum on each plate or test strip series.
12. Use only distilled deionised or any similar high quality water for preparation of reagents.
13. When preparing the buffers, etc., measure the required volume.
14. The Stop Solution contains sulphuric acid, which is corrosive.\*
15. All unused biological materials should be disposed according to the local, regional and national regulations.

## Recommendations!

There is always a surplus volume of liquid reagents. The volume mentioned on the label is the minimum obtainable.

Strips with broken seal can be stored at 2-8°C (36-46°F) for up to 4 weeks.

## Procedure

1. All reagents should equilibrate to room temperature 18-25°C (64-77°F) before use. Label each strip with a number.
2. Add samples  
The provided negative and positive control sera are used for both serum and milk testing.

### Serum Samples

- A. Add 90 µL of Sample Dilution Buffer to each well that will be used for serum samples and serum controls.
- B. Add 10 µL of Positive Control Serum (Reagent A) and 10 µL of Negative Control Serum (Reagent B) respectively to selected wells coated with BVDV viral antigen. For confirmation purposes it is recommended to run the control sera in duplicates.
- C. Add 10 µL of serum sample to a selected well coated with BVDV viral antigen. The samples can be tested in singlicates or in duplicates. However, for confirmation purposes it is recommended to run the samples in duplicates.

Continue at step #3.

### Milk Samples

- A. For addition of controls, see "Serum Samples" (point A and B).
- B. Add 100 µL of skim milk sample to a selected well coated with BVDV viral antigen. The samples can be tested in singlicates or in duplicates. However, for confirmation purposes it is recommended to run the samples in duplicates.

Continue at step #3.

3. Shake the plate thoroughly on a plate shaker for a few seconds. Seal the plate/strips and incubate at 37°C (98.6°F) for 1 hour, or at 2-8°C (36-46°F) over night (16-20 hours).
4. Rinse the plate/strips 3 times with PBS-Tween Buffer. At each rinse cycle fill up the wells, empty the plate and tap hard to remove all remains of fluid.
5. Add 100 µL of HRP Conjugate to each well. Seal the plate/strips and incubate at 37°C (98.6°F) for 1 hour.
6. Repeat step #4.
7. Add 100 µL Substrate Solution to each well. Incubate for 10 minutes at room temperature 18-25°C (64 -77°F). Begin timing when the first well is filled.
8. Stop the reaction by adding 50 µL of Stop Solution to each well and mix thoroughly. Add the Stop Solution in the same order as the Substrate Solution in step #7.
9. Measure the optical density (OD) of the controls and samples at 450 nm in a microplate photometer (use air as blank). Measure the OD within 15 minutes after the addition of Stop Solution to prevent fluctuation in OD values.

## Calculations

Calculations of results are done as described below.

### Calculation of Percent Positivity Values (PP)

All OD Values for the test samples as well as the Negative Control are related to the OD value of the positive control as follows:

$$PP = \frac{OD_{\text{Sample or Negative control}}}{OD_{\text{Positive control}}} \times 100$$

## Interpretation of the results

### Criteria for test validity

To ensure validity, the duplicate OD values of the positive control should not differ more than 25 % from the mean value of the two duplicates.

Additionally, the control values should fall within the following limits:

**OD Positive Control** > 0.7

**PP Negative Control, short incubation** < 17

**PP Negative Control, long incubation** < 30

Should any of these criteria not be fulfilled, the test is invalid. For invalid tests, technique may be suspected and the assay should be repeated.

If samples are tested in duplicates, each separate duplicate PP value should be equally interpreted, *i.e.* positive or negative. In case of discrepancy it is recommended to retest the sample.

It is recommended to re-test all samples with a positive reaction with an ELISA with control antigen line, to eliminate any non-specific reacting samples *i.e.* the SVANOVIR® BVDV-Ab ELISA, Art. No. SV-104881 and SV-104882.

## Interpretation

### Short incubation (1 hour incubation)

Sample	PP	Interpretation
Serum	< 20 ≥ 20	Negative Positive
Milk	< 12 ≥ 12	Negative Positive

### Long incubation (over night incubation)

Sample	PP	Interpretation
Serum	< 30 ≥ 30	Negative Positive
Milk	< 17 ≥ 17	Negative Positive

## References

1. Barber, D.M.L., Nettleton, P.F., and Herring, J.A. (1985) Disease in a dairy herd associated with the introduction and spread of bovine virus diarrhoea virus. *Vet. Rec.* 2, 459-464.
2. Van Oirschot, J.T. (1983) Congenital infections with nonarbo togaviruses. *Vet. Microbiol.* 8, 321-361.
3. Howard, C.J., Brownlie, J. and Clarke, M.C. (1987). Comparison by the neutralization assays of pairs of noncytopathogenic and cytopathogenic strains of bovine virus diarrhoea virus isolated from cases of mucosal disease. *Vet. Microbiol.* 13, 361-369.
4. Reggiardo, C., and Kaeberle, M.L. (1981) Detection of bacteremia in cattle inoculated with bovine viral diarrhoea virus. *The Am. J. of Vet. Res.* 42, 218-221.

### **\*DANGER: Stop solution (sulphuric acid)**

May be corrosive to metals. Causes skin irritation. Causes serious eye irritation.

Keep only in original container. Wear eye protection/ face protection. Wear protective gloves.

**IN CASE OF CONTACT WITH EYES:** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/ physician. If eye irritation persists: Get medical advice/ attention.

**IN CASE OF CONTACT WITH SKIN:** Wash with plenty of soap and water. Take off contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention. Absorb spillage to prevent material damage.

INDICAL BIOSCIENCE GmbH  
Deutscher Platz 5b  
04103 Leipzig  
Germany  
[www.indical.com](http://www.indical.com)

**Customer Service**  
[support@indical.com](mailto:support@indical.com)

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## Referenzen

1. Barber, D.M.L., Nettleton, P.F., and Herring, J.A. (1985) Disease in a dairy herd associated with the introduction and spread of bovine virus diarrhoea virus. Vet. Rec. 2, 459-464.
2. Van Oirschot, J.T. (1983) Congenital infections with nonarbo togaviruses. Vet. Microbiol. 8, 321-361.
3. Howard, C.J., Brownlie, J. and Clarke, M.C. (1987). Comparison by the neutralization assays of pairs of noncytopathogenic and cytopathogenic strains of bovine virus diarrhoea virus isolated from cases of mucosal disease. Vet. Microbiol. 13, 361-369.
4. Reggiardo, C., and Kaeberle, M.L. (1981) Detection of bacteremia in cattle inoculated with bovine viral diarrhoea virus. The Am. J. of Vet. Res. 42, 218-221.

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