SVANOVIR® PRV gE-Ab

Pseudorabies/Aujeszky's Disease
Antibody Test

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Contents	Art. No. SV-104902	Art. No. SV-104904
Microtitre plate Microtitre plates (96 wells) coated with non-infectious PRV antigen (sealed and stored dry)	2 (Strips) 12 x 8	10 (Plates)
Conjugate Ready-to-use (horseradish peroxidase conjugated anti-PRV gE monoclonal antibodies)	1 x 24 mL	1 x 120 mL
PBS-Tween Solution 20 x concentrate	1 x 125 mL	3 x 125 mL
Sample Dilution Buffer	1 x 8 mL	1 x 40 mL
Substrate Solution (Tetramethylbenzidine in substrate buffer containing H ₂ O ₂) - STORE IN THE DARK!	1 x 20 mL	1 x 100 mL
Stop Solution Contains sulphuric acid (2M) - DANGER!	 1 x 10 mL	2 x 25 mL
A. Positive Control Serum - Contains preservatives	1 x 2.5 mL	1 x 4.5 mL
B. Negative Control Serum - Contains preservatives	1 x 2.5 mL	1 x 4.5 mL

This manual covers the following SVANOVIR® PRV gE-Ab ELISA kits:
Article number SV-104902 and SV-104904

Pseudorabies/Aujeszky's Disease

Antibody Test

Name and Application

SVANOVIR® PRV gE-Ab is an Enzyme Linked Immunosorbent Assay (ELISA) for the detection of PRV gE antibodies in porcine serum and plasma samples. The test will discriminate between animals infected and animals vaccinated with a gE-deleted vaccine).

General information

Pseudorabies Virus (PRV) or suid herpesvirus¹, first isolated by Aujeszky (1902) from an ox with signs of rabies, has the pig as the primary host, while cattle and some other animal species are secondary hosts with infected pigs as the initial source of infection. In countries with a high prevalence of endemic pseudorabies, this disease, being often, particularly in piglets, fatal, causes severe losses annually¹. In piglets, but also in older pigs, neurological disorders are intermingled with symptoms from the respiratory tract. The disease is less pronounced in older pigs and, after recovery, *i.e.* in adult pigs, a latent infection persists lifelong. From such asymptomatic pigs, PRV has been isolated from trigeminal ganglia and tonsils². Virus transmission occurs mainly through nasal secretions by contact or even by the aerogenic route¹. The infection of pregnant sows may reach via the placenta the fetuses and cause abortion, mummification or stillbirth. All infected herds in endemic regions should be monitored for the presence of infection and uninfected herds protected by control measures. Some countries practice vaccination, while some others try to control the spread by culling seropositive pigs or by combination of both measures.

Differentiation of vaccinated from infected pigs

If both control measures, *i.e.* vaccination and eradication by culling seropositive pigs are (or are to be) introduced, a differentiation of vaccinated pigs from those infected must be introduced, as both infected and vaccinated pigs are seropositive. Most vaccinal strains used are well-known attenuated strains with a deletion in their genome (DNA) of the gene sequence coding for the virulence factor, the glycoprotein E, (gE/gI)². Consequently, vaccinated pigs lack antibodies against the gE antigen. The differentiation is thus based on the presence of antibodies to gE in infected pigs but the absence of these antibodies in vaccinated pigs³.

Since both the vaccinal and field strains of PRV possess another glycoprotein, gB (gII), antibodies to gB are present in sera of both infected and vaccinated pigs but absent in sera of non-infected and non-vaccinated pigs⁴. The SVANOVIR® PRV gE and gB ELISA systems are based on these premises.

Principle

The kit procedure is based on the blocking Enzyme Linked Immunosorbent Assay (ELISA). In this procedure the gE-test, microtitre plates are coated with PRV antigen. This antigen is derived from an inactivated PRV strain, known to possess the glycoprotein gE (as well as gB). The indicator used for the detection of presence or absence of antibodies against gE is a monoclonal antibody (mAb) specifically directed against gE. The mAb are conjugated with the enzyme horse-radish peroxidase (HRP). If antibodies are absent in the test serum, *i.e.* from non-infected, non-vaccinated or vaccinated pigs, the gE antigen will remain free and the conjugated mAb subsequently added will bind. The enzyme present (HRP) will catalyze the

added substrate, which results in a colour reaction. If, on the other hand, antibodies against gE are present in the serum, *i.e.* from pigs infected with field virus, these antibodies will bind to the gE antigen in the well. This will block the possibility of the conjugated mAb subsequently added to react with the gE antigen. The unbound conjugated mAb will be removed by the washing procedure followed, resulting in no colour reaction when the substrate is added.

Materials needed but not provided

1. Precision pipettes
2. Disposable pipette tips
3. Distilled, deionised or any similar high quality water
4. Wash bottle, multichannel pipettor or plate washer
5. Container: 1 to 2 litres for PBS-Tween
6. Microplate photometer, 450 nm filter

Specimen information

100 µL of blood serum or plasma is required for each sample well. Fresh, refrigerated or previously frozen serum or plasma may be tested.

Preparation of reagents

PBS-Tween Buffer:

Dilute the PBS-Tween Solution 20 x concentrate 1/20 in distilled water. Prepare 500 mL per plate by adding 25 ml PBS-Tween solution to 475 ml distilled water and mix thoroughly.

N.B. Please check that there is no crystal precipitation in the bottle. If crystals are seen, please warm and shake well.

Precautions

1. Carefully read and follow all instructions.
2. Store the kit and all reagents at 2-8°C (36-46°F).
3. All reagents should equilibrate to room temperature 18-25°C (64-77°F) before use.
4. Handle all materials according to the Good Laboratory Practice.
5. Do not mix components or instruction manuals from different test kit batches.
6. Care should be taken to prevent contamination of kit components.
7. Do not use test kit beyond date of expiry.
8. Do not eat, drink, or smoke where specimens or kit reagents are handled.
9. Use a separate pipet tip for each sample.
10. Do not pipet by mouth.
11. Include positive and negative serum on each plate or test strip series.
12. Use only distilled deionised or any similar high quality water for preparation of reagents.
13. When preparing the buffers, etc., measure the required volume.
14. The Stop Solution contains sulphuric acid, which is corrosive.*
15. All unused biological materials should be disposed according to the local, regional and national regulations.

Recommendations!

For article SV-104902, strips with broken seal can be stored at 2-8°C (36-46°F) for up to 4 weeks.

There is always a surplus volume for the liquid reagent.

The volume mentioned on the label is the minimum obtainable.

Procedure

1. All reagents should equilibrate to room temperature 18-25°C (64-77°F) before use. Label each strip with a number.

Short sample incubation:

2. Add 100 µL of Positive Control Serum (reagent A) and 100 µL of Negative Control Serum (reagent B) respectively, to selected wells coated with PRV antigen. For confirmation purposes it is recommended to run the control sera in duplicates.
3. Add 100 µL of undiluted serum sample to a selected well coated with PRV antigen. For confirmation purposes it is recommended to run the samples in duplicates.
4. Seal the plate/strip and incubate at 37°C (98.6°F) for 1 hour.
Subsequently continue as from step #8.

Over Night Incubation:

5. Add 25 µL of Sample Dilution Buffer to all wells.
6. Repeat step #2 and step #3.
7. Seal the plate/strip and incubate at room temperature, 18-25°C (64-77°F) over night (12-18 hours).
8. Rinse the plates/strips 3 times with PBS-Tween Buffer: at each rinse cycle fill the wells, empty the plate and tap hard to remove all remains of fluid
9. Add 100 µL of Conjugate to each well. Seal the plate/strips and incubate at 37°C (98.6°F) for 30 minutes or at room temperature 18-25°C (64-77°F) for 60 minutes.
10. Repeat step #8.
11. Add 100 µL of Substrate Solution to each well. Incubate for 15 minutes at room temperature 18-25°C (64-77°F). Begin timing when the first well is filled
12. Stop the reaction by adding 50 µL of Stop Solution to each well and mix thoroughly. Add the Stop Solution in the same order as the Substrate Solution in step #11.

13. Measure the optical density (OD) of the controls and samples at 450 nm in a microplate photometer (use air as a blank). Measure the OD within 15 minutes after the addition of Stop Solution to prevent fluctuation in OD values.

Calculations

1. Calculate the mean optical density (OD) value for each of the control sera.
2. Calculate the optical density percentage (ODp) for each sample and for the positive control sera using the following formula:

$$\text{ODp} = \frac{\text{OD}_{\text{Sample/Pos control}}}{\text{OD}_{\text{Negative control}}} \times 100$$

Interpretation of the results

Criteria for test validity

To ensure validity, the duplicate OD values of the negative control should not differ more than 25% from the mean value of the two duplicates. Additionally, the control values should fall within the following limits:

OD Negative Control 0.8 -1.6

ODp Positive Control < 30

Should any of these criteria not be fulfilled, the test is invalid. For invalid tests, technique may be suspect and the assay should be repeated.

Interpretation

Short sample incubation:

1. **ODp \leq 50 = positive**

Serological diagnosis:
Field virus infected.

2. **ODp $>$ 60 = negative**

Serological diagnosis: Either

a) non-infected with field virus and no vaccinated against Aujeszky or

b) non-infected with field virus but vaccinated with an Aujeszky disease gE-deleted vaccine.

3. **ODp 51 - 60 = doubtful**

Retest the sample. If the test result is still doubtful, a new sample should be taken from the animal and both samples should be run in parallel.

Over Night sample incubation:

1. **ODp \leq 45 = positive**

Serological diagnosis:
Field virus infected.

2. **ODp $>$ 55 = negative**

Serological diagnosis: Either

a) non-infected with field virus and no vaccinated against Aujeszky or

b) non-infected with field virus but vaccinated with an Aujeszky disease gE-deleted vaccine.

3. **ODp 46 - 55 = doubtful**

Retest the sample. If the test result is still doubtful, a new sample should be taken from the animal and both samples should be run in parallel.

To evaluate whether the animal, for which a test serum score negative, is vaccinated, it is recommended to test the serum in the SVANOVIR® PRV gB-Ab test (Art. No. SV-104905/SV-104906).

References

1. Gustafson, D.P. (1986) Pseudorabies. Diseases of Swine 6th Edition. Edited by A.D. Leman, R.D. Glock, W.C. Mengeling, R.H. Penny, E. Scholl, and B. Straw. Iowa State Uni. Press, pp. 274-289.
2. Rziha, H.-J., Döller, P.C., and Wittman, G. (1982) Detection of Aujeszky's Disease virus and viral DNA in tissues of latently infected pigs. Current Topics in Veterinary Medicine and Animal Science., Vol. 17, Aujeszky's Disease. Edited by G. Wittman and S.A. Hall. Martinus Nijhoff, The Hague, pp 205-214.
3. Baskerville A., Mc Ferran, J.B., and Dow, C. (1973) Aujeszky's Disease in pigs. The Vet. Bull. 43,465-480.
4. Soerensen, K.J., and Lei, J.C. (1986) Aujeszky's Disease: blocking ELISA for the detection of serum antibodies. J. Vir. Mth. 13, 171-181.

***DANGER: Stop solution (sulphuric acid)**

May be corrosive to metals. Causes skin irritation. Causes serious eye irritation.

Keep only in original container. Wear eye protection/ face protection. Wear protective gloves.

IN CASE OF CONTACT WITH EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/ physician. If eye irritation persists: Get medical advice/ attention.

IN CASE OF CONTACT WITH SKIN: Wash with plenty of soap and water. Take off contaminated clothing and wash it before reuse. If skin irritation occurs: Get medical advice/attention. Absorb spillage to prevent material damage.

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Customer Service
support@indical.com

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Referencias

1. Gustafson, D.P. (1986) Pseudorabies. Diseases of Swine 6th Edition. Edited by A.D. Leman, R.D. Glock, W.C. Mengeling, R.H. Penny, E. Scholl, and B. Straw. Iowa State Uni. Press, pp. 274-289.
2. Rziha, H.-J., Döller, P.C., and Wittman, G. (1982) Detection of Aujeszky's Disease virus and viral DNA in tissues of latently infected pigs. In Current Topics in Veterinary Medicine and Animal Science., Vol. 17, Aujeszky's Disease. Edited by G. Wittman and S.A. Hall. Martinus Nijhoff, The Hague, pp 205-214.
3. Baskerville A., Mc Ferran, J.B., and Dow, C. (1973) Aujeszky's Disease in pigs. The Vet. Bull. 43,465-480.
4. Soerensen, K.J., and Lei, J.C. (1986) Aujeszky's Disease: blocking ELISA for the detection of serum antibodies. J. Vir. Mth. 13, 171-181.

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INDICAL BIOSCIENCE GmbH