

IndiMag[®] 2

Installation Guide

The IndiMag 2 requires a clean level surface with minimum dimensions of 65 cm in width, 75 cm in height and 75 cm in depth. The installation site should be free of vibrations, strong magnetic fields, direct sunlight or UV light, drafts, excessive moisture, excessive dust, and large temperature fluctuations.

Maintain a minimum clearance of 5 cm to each side of the IndiMag 2. Do not block any airflow beneath the unit. Slits and openings that ensure the instrument ventilation must not be covered.

Further information and support

- IndiMag 2 Handbook: www.indical.com/indimag2
- Technical assistance: support@indical.com

Unpacking the instrument

1. Remove the lid from the box (see Figure 1).
2. Remove the accessories box.
3. Remove the upper packing foam and plastic bag.
4. Remove the cardboard sleeve from the instrument.

Note: It is recommended to retain the packaging.

5. Lift the instrument and place it on a benchtop.

Note: The instrument weights 32 kg and requires at least two people to lift.



Figure 1. Unpacking the IndiMag 2 instrument and lifting the device by two people.

First instrument setup

1. Remove the transport locks (see Figure 2a and b). Press the button on the locking pins and pull out simultaneously.
Note: Store the transport locks in the designated spaces inside the instrument (see Figure 2c). Press the button on the locking pins and simultaneously push into designated slots.
2. Adjust the four feet to level the instrument on your installation site (see Figure 2d).
3. Connect the mains cable.
4. Turn on the IndiMag 2.




Figure 2. Positions of transport locks and instrument feet.

For up-to-date licensing information and product-specific disclaimers, see the respective INDICAL kit handbook or user manual. Trademarks: IndiMag® and INDICAL® (INDICAL BIOSCIENCE). Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are not to be considered unprotected by law.
HB-2660-EN-001 04/2024 © 2024 INDICAL BIOSCIENCE GmbH, all rights reserved.