

# virotype<sup>®</sup> CaPV pan Positive Control

May 2024 | EN

**REF** 30 reactions (cat. no. PC286745)



INDICAL BIOSCIENCE GmbH, Deutscher Platz 5b, 04103 Leipzig, Germany

**virotype CaPV pan Positive Control**

**(30)**

**Cat. no.**

**PC286745**

**Number of reactions**

**30**

CaPV pan Positive Control

1 x 240 µl

## Intended use

The virotype CaPV pan Positive Control contains synthetically derived DNA from *Capripoxvirus* (CaPV). This product can be used as a positive control in real-time PCR.

**For research use only. This product is not intended for the diagnosis, prevention or treatment of a disease.**

## Preparation of PCR setup

Before use, thaw the reagent. Invert the tube several times to mix, and spin down briefly.

It is recommended to use 8 µl of the virotype CaPV pan Positive Control per 25 µl PCR reaction.

## Storage and handling

Upon receipt, store the material at –30°C to –15°C, away from any sources of contaminating DNA or RNA, especially amplified DNA products and material used for PCR master mixes. Protect from light during storage and handling. Use aerosol barrier pipet tips for pipetting. Do not use after the expiration date printed on the label. Avoid repeated thawing and freezing because this may reduce

performance. If used intermittently, aliquot the material and freeze.

## Quality Control

Each lot of this product was manufactured and released in accordance with INDICAL's ISO-certified Quality Management System.

## Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles.

For more information, please consult the appropriate safety data sheets (SDSs). These are available from your regional sales manager or at [compliance@indical.com](mailto:compliance@indical.com).

All sample residues and objects that have come into contact with samples must be decontaminated or disposed of as potentially infectious material.

## Warning and precautions

Dispose of all materials in accordance with national and local regulations.

## Change index

Product Sheet	Version	Change
HB-2656-EN-001	May 2024	Product implementation

**Trademarks:** virotype<sup>®</sup> (INDICAL BIOSCIENCE GmbH). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

HB-2656-EN-001 © 2024 INDICAL BIOSCIENCE GmbH, all rights reserved.

**INDICAL**  
BIOSCIENCE

Ordering: [www.indical.com/contact](http://www.indical.com/contact)  
Technical Support: [support@indical.com](mailto:support@indical.com)  
Website: [www.indical.com](http://www.indical.com)