

APPLICATION NOTE

Rapid detection using the IndiMag® Pathogen Kit and the virotype® IHNV/VHSV RT-PCR Kit

Daland C. Herrmann, Marco Labitzke, Douglas Marthaler, Oliver Sasse, Stefanie Fritsche, and Ansgar Flammersfeld

INDICAL BIOSCIENCE GmbH, Deutscher Platz 5b, 04103 Leipzig, Germany



Introduction

Infectious Hematopoietic Necrosis Virus (IHNV) and **Viral Hemorrhagic Septicemia Virus (VHSV)** are highly contagious RNA viruses that spread systemically and cause hematopoietic necrosis and hemorrhagic lesions, leading to high mortality rates and economic losses in salmonid species, particularly affecting rainbow trout and other farmed fish [1].

Effective monitoring and rapid diagnosis are critical to prevent outbreaks and ensure biosecurity [2]. Real-time Reverse Transcription Polymerase Chain Reaction (real-time RT-PCR), provides a fast, sensitive, and highly specific method for detecting viral RNA in fish tissues and is key to controlling the spread of these viruses, maintaining healthy fish populations, and protecting the aquaculture industry [2, 3].

In this application note, INDICAL's **IndiMag Pathogen Kit** for DNA/RNA extraction, coupled with the **virotype IHNV/VHSV RT-PCR Kit** illustrate a reliable and efficient solution for the molecular detection of *Infectious Hematopoietic Necrosis Virus* and *Viral Hemorrhagic Septicemia Virus* to support disease surveillance and management in fish farming.

Keywords: Biosecurity, pathogen detection, viral RNA, IHNV (*Infectious Hematopoietic Necrosis Virus*), VHSV (*Viral Hemorrhagic Septicemia Virus*), RNA extraction, real-time RT-PCR

Materials and Methods

Sample preparation and RNA extraction

IHNV-positive (VHSV-negative) fish samples from a recent outbreak were kindly provided by a Veterinary State Laboratory in Germany. Organs (brain, kidney, liver, and spleen) were dissected, minced into small pieces and pooled into organ homogenates.

The organ homogenates were divided to enable a direct comparison between the **IndiMag Pathogen Kit** (INDICAL BIOSCIENCE) with **Pretreatment T1** and the **Reference Method** according to the official protocol published in the "Amtliche Methodensammlung" (Official collection of analytical methods) by the Friedrich-Loeffler-Institut (FLI) using silica spin columns for RNA extraction [4, 5, 6].

The protocol with the IndiMag Pathogen Kit with Pretreatment T1 was as follows:

1. 25 mg of the organ mix was transferred into 2 mL tubes containing 2.38 mm metal PowerBeads® (QIAGEN) and 300 µL of Phosphate-buffered saline (PBS) was added.
2. The samples were homogenized for 1 minute at 6.5 m/s using the FastPrep®-24 instrument (MP Biomedicals, LLC) and centrifuged for 2 minutes at 4,000 × g.
3. Then, 200 µL of the supernatant was used for RNA extraction with the IndiMag Pathogen Kit on the KingFisher™ Flex (Thermo Fisher Scientific Inc.).

The protocol using the Reference Method was as follows:

1. The organ mix was homogenized using a mortar and pestle, transferred into a 2 mL tube, and centrifuged for 15 minutes at 2,000 × g.
2. Then, 70 µL of the supernatant was mixed with 70 µL PBS and 560 µL Buffer AVL (QIAGEN), briefly vortexed, and incubated at room temperature for 15 minutes.
3. The RNA was extracted by silica column according to the manufacturer's protocol.

The IndiMag Pathogen Kit is available in two formats:

- Manual-fill where the user adds all reagents to the required plates.
- Prefilled format for the IndiMag 2, IndiMag 48s, and KingFisher™ systems, where reagents are pre-dispensed into plastics and only the sample and Lysis Buffer need to be added.

The IndiMag Pathogen Kit facilitates lysis, binding, and purification of total nucleic acids.

Detection and amplification

The extracted nucleic acids were analyzed with the **virotype IHN/VHSV RT-PCR Kit** (INDICAL BIOSCIENCE), following the manufacturer's instructions [7]. This ready-to-use triplex RT-qPCR assay enables the simultaneous identification and differentiation of RNA from IHN (Cy5) and VHSV (FAM) across all major European, North American, and Pacific genotypes, and includes an endogenous salmonid β -actin control (HEX) to support result interpretation.



Figure 1: virotype IHN/VHSV RT-PCR Kit (Cat. no. VT290005)

The virotype IHN/VHSV RT-PCR Kit is approved by the Friedrich-Loeffler-Institut (FLI-C 133), Germany's national reference laboratory for animal health and zoonoses.

In brief, the PCR reactions (25 µL total volume, including 5 µL RNA input and 20 µL of Master Mix) were performed on a Bio-Rad® CFX96™ PCR cycler (Bio-Rad, Hercules, CA, USA) [7].

The analytical sensitivity of the virotype IHN/VHSV RT-PCR Kit was previously validated, with a limit of detection of 2 copies per reaction for IHN and VHSV.

The correlation coefficient and PCR efficiency are 0.987 and 110.1 % for IHN and 0.994 and 111.2 % for VHSV, respectively [8].



Representative application results

The virotype IHN/VHSV RT-PCR Kit was employed to compare the extraction chemistries of the IndiMag Pathogen Kit with Pretreatment T1 and the Reference Method, using IHN-positive (VHSV-negative) fish tissue homogenate samples.

In this dataset, the IndiMag Pathogen Kit and the Reference Method equivalently extracted genomic nucleic acids, with both methods illustrating strong and consistent amplification of the salmonid β -actin EC indicating efficient and reproducible RNA recovery. The IndiMag Pathogen Kit with Pretreatment T1 yielded a lower mean C_T of 2.27 for IHN-RNA identification, suggesting improved extraction performance compared to the Reference Method (Table 1). The samples tested negative for VHSV with the virotype IHN/VHSV RT-PCR Kit.

Table 1: Comparative C_T values of IHN-RNA and the endogenous Internal Control (EC) extracted using the IndiMag Pathogen Kit with Pretreatment T1 vs. the Reference Method. C_T values are shown for IHN-positive fish organ homogenate samples. All samples tested negative for VHSV (data not shown).

Material and sample number	Reference Method		IndiMag Pathogen Kit	
	IHN C_T	EC C_T	IHN C_T	EC C_T
Fish organ homogenate #1	33.80	19.17	33.73	19.05
Fish organ homogenate #2	24.32	16.73	21.10	18.16
Fish organ homogenate #3	18.00	16.68	15.32	18.33
Fish organ homogenate #4	20.09	17.70	18.87	17.26
Fish organ homogenate #5	22.40	16.94	19.97	17.86
Fish organ homogenate #6	38.42	16.25	32.80	17.07
Fish organ homogenate #7	19.16	18.92	14.92	17.63
Fish organ homogenate #8	16.47	15.30	15.51	16.37
Fish organ homogenate #9	18.65	15.07	14.70	16.29
Fish organ homogenate #10	17.47	16.06	16.15	15.93
Fish organ homogenate #11	16.87	15.87	15.40	15.81
Fish organ homogenate #12	26.71	16.06	25.58	16.65
Fish organ homogenate #13	27.62	16.53	26.53	16.26
Mean C_T	23.08	16.75	20.81	17.13




Figure 2: Graphical comparison of C_T values of the IHN-positive fish organ homogenate samples 1-13 obtained using the IndiMag Pathogen Kit method with Pretreatment T1 and Reference Method as shown in Table 1.

Conclusion

INDICAL's IndiMag Pathogen Kit with Pretreatment T1 workflow delivers consistent, high-quality nucleic acid extraction from fish tissue samples, ensuring robust performance comparable to the Reference Method. Specifically designed for tissue matrices, it releases more viral material into the supernatant, providing a reliable alternative to the gold standard. In addition to the manual-fill format, the IndiMag Pathogen Kit is also available in prefilled versions for automated nucleic acid purification platforms such as the IndiMag 2, IndiMag 48s, KingFisher, and equivalent systems, further enabling efficient processing of large sample numbers with reduced hands-on time and increased laboratory capacity.

The virotype IHN/VHSV RT-PCR Kit provides a highly sensitive and specific tool for the reliable detection of IHN and VHSV in aquaculture samples. It covers all major genotypes (European, North American, and Pacific) and demonstrates no cross-reactivity with other aquatic viruses. The assay is officially approved by the FLI (FLI-C 133), Germany's Federal Research Institute for Animal Health [4, 7].

Comprehensive validation data [8] are available upon request to support regulatory or operational implementation.

References

1. Hoferer M, Akimkin V, Skrypski J, Schutze H, Sting R. Improvement of a diagnostic procedure in surveillance of the listed fish diseases IHN and VHS. J Fish Dis. 2019;42(4):559–72.
2. WOAH (World Organisation for Animal Health). Manual of Diagnostic Tests for Aquatic Animals. 24th ed., Paris, France. 2024. Chapter 2.3.5: Infection with infectious haematopoietic necrosis; Chapter 2.3.10: Infection with viral haemorrhagic septicaemia virus.
Available online at: <https://www.woah.org/en/what-we-do/standards/codes-and-manuals/aquatic-manual-online-access/>
3. European Commission. Commission Implementing Decision 2015/1554 of 11 September 2015 laying down rules for the application of Directive 2006/88/EC as regards requirements for surveillance and diagnostic methods. Official Journal of the European Union. 2015;L 247:1–62.
4. Friedrich-Loeffler-Institut (FLI). Infektiöse Hämatopoetische Nekrose (IHN) und Virale Hämorrhagische Septikämie (VHS) der Salmoniden: Amtliche Methode und Falldefinition. Amtliche Methodensammlung und Falldefinitionen: Anzeigepflichtige Tierseuchen. Greifswald – Insel Riems, Germany. 2025
Available online at: <https://www.fli.de/de/publikationen/amtliche-methodensammlung>
5. INDICAL BIOSCIENCE (INDICAL). IndiMag Pathogen pretreatment protocol for mechanical disruption of tissue samples.
Internal document, available upon request from manufacturer.
6. INDICAL BIOSCIENCE (INDICAL). IndiMag Pathogen Kit Handbook.
Available online at: shop.indical.com
7. INDICAL BIOSCIENCE (INDICAL). virotype IHN/VHSV RT-PCR Kit Handbook.
Available online at: shop.indical.com
8. INDICAL BIOSCIENCE (INDICAL). virotype IHN/VHSV RT-PCR Kit Validation Report.
Internal document, available upon request from manufacturer.

Access these and other related INDICAL protocols on the respective INDICAL BIOSCIENCE product pages. For a complete overview of INDICAL's solutions for nucleic acid extraction, detection, and instrumentation, please visit shop.indical.com or contact your local INDICAL representative.

Website: www.indical.com | Email: support@indical.com | Call: +49 341 124 54 0
Email (Americas): us_support@indical.com | Call (Americas): +1 (407) 766 5128
Connect with us on LinkedIn: www.linkedin.com/company/indical

Disclaimer: This application note provides an illustrative example of using the IndiMag Pathogen Kit for DNA/RNA purification, including pretreatment methods and related INDICAL solutions such as virotype assays. It highlights defined protocols that demonstrate the versatility of the reagents and workflows. However, this document does not constitute a fully validated protocol. Each laboratory is responsible for performing its own validation in accordance with applicable regulations and institutional guidelines. Only the specifications and information provided in the official user manuals apply. For up-to-date licensing information and product-specific disclaimers, see the relevant kit handbook or user manual. Assays for veterinary use only. Reagents for research use only, not for use in diagnostic procedures. INDICAL BIOSCIENCE makes no warranties or representations regarding the accuracy, reliability, or suitability of the application for any particular purpose and assumes no liability for any use or consequences arising from this information, including errors or omissions. All other terms and conditions remain unaffected. Regulatory requirements vary by country; products may not be available in your geographic area. Product images may differ from the actual product. This document is not intended to encourage the use of any products in ways that might infringe upon the intellectual-property rights of others.

Trademarks: INDICAL®, IndiMag®, and virotype® (INDICAL BIOSCIENCE GmbH); Bio-Rad® and CFX96™ (Bio-Rad Laboratories, Inc.); Cy5® (GE HealthCare); FastPrep® (MP Biomedicals, LLC); FAM™, HEX™, and KingFisher™ (Thermo Fisher Scientific Inc.); PowerBead®, and QIAGEN® (QIAGEN Group). Registered names, trademarks, etc., even if not explicitly marked as such, may still be protected by law. All other trademarks are the property of their respective owners. The use of brand or product names is solely for identification purposes and does not imply endorsement.

IN-AppNote-2501-EN-Oct-2025 © 2025 INDICAL BIOSCIENCE, all rights reserved.