

Check&Trace **Salmonella 2.0** QS5

For use with the ThermoFisher QuantStudio™ 5 Real-Time PCR System (96-well 0.1 mL or 0.2 mL Block)



1 INTENDED USE

The Check&Trace Salmonella 2.0 (CTS 2.0) assay is a qualitative, semi-automated real-time PCR test designed for the confirmation and typing of presumptive *Salmonella* isolates from enriched culture media. The confirmation procedure advances a suspected (presumptive) *Salmonella* result (according to ISO-6579-1) to a confirmed *Salmonella spp.* result or a *Salmonella* negative result. The typing procedure generates a serovar result or a “genovar” code both uniquely defining each bacterial isolate tested.

Testing is performed after a simple sample preparation step of bacteria from culture. The sample is then used as input material for a 6-tube multiplex real-time PCR assay performed on the ThermoFisher QuantStudio™ 5 Real-Time PCR System. Presence and absence of target sequences for confirmation and typing are assessed for each multiplex PCR and a final test result is generated automatically using dedicated software.

CTS 2.0 is intended for (as mentioned in ISO 6579):

- Products intended for human consumption and the feeding of animals,
- Environmental samples in the area of food production and food handling,
- Samples from the primary production stage such as animal faeces, dust, and swabs.

2 SUMMARY AND EXPLANATION OF THE PROCEDURE

Salmonella is one of the main causes of food poisoning. Detection of *Salmonella* in food and feed stuff is therefore mandatory in many countries. Standard protocols for detection of *Salmonella* are time-consuming, taking several days to generate a final positive or negative test result. (ISO-6579/1). If positive for *Salmonella*, many countries require further typing of *Salmonella* for epidemiological and health risk assessments. Traditionally this is done using serotyping for the somatic (O) lipopolysaccharide and flagellar (H) protein antigens (ISO-6579/3). However, *Salmonellae* are antigenically complex with over 2 400 serovars having various combinations of the 46 O-antigens and 85 H-antigens (3). So, once positive, further typing of a positive *Salmonella* may again take several days to complete requiring many O- and H-antisera for full typing.

DNA typing differs from serotyping. With serotyping the presence of antigens on the cell surface and flagella are detected. This is based on expression of genes located on two specific segments of the *Salmonella* genome. The CTS 2.0 assay detects genetic variation at 21 loci scattered over the whole *Salmonella* genome: this generates specific *Salmonella* genotypes, also called Genovars. The CTS 2.0 database links Genovars to a collection of well-characterized *Salmonella* Serovars. A Genovar will be given as test result if there is no established Serovar for the Genovar generated by the CTS 2.0 assay.

With the CTS 2.0 assay confirmation and typing of *Salmonella* can be performed in approximately 2.5 hours starting from bacteria from a XLD or Nutrient Agar (NA) plate. The presence of DNA target sequences for typing and confirmation are detected in parallel with various internal controls. CTS 2.0 uses a quick and simple sample preparation process and minimizes operator intervention once the samples are placed onto the real-time PCR system.

PRINCIPLES OF THE PROCEDURE

A worklist is generated on the CTS 2.0 analysis portal. Bacterial specimens are collected from XLD or NA media using “colony samplers” included in the CTS 2.0 kit. Colony samplers are then placed in 1.5 mL tubes prefilled with 200 µL Sample Buffer. Bacteria are resuspended in this buffer and tubes are closed and transferred to a heating block pre-heated at 98 °C. Tubes are incubated for 10 minutes at 98 °C and allowed to cool down to room temperature for 5 minutes. This procedure generates a crude lysate

from the specimens releasing the bacterial DNA in the Sample Buffer. This crude lysate is suitable for the CTS 2.0 real-time amplification process without further purification.

140 µL of the crude lysate is transferred to a fresh tube prefilled with 35 µL PCR Master Mix and mixed by pipetting up and down 3 times. Next 25 µL of the crude lysate/Master Mix solution is added to each of 6 tubes of the CTS 2.0 PCR Reagent Strip. This Reagent Strip contains all PCR primers and probes for target sequence amplification and detection using real-time PCR. Finally, the PCR Reagent Strip is placed in the real-time PCR instrument and the amplification process is started.

The amplified DNA targets are detected using hydrolysis (TaqMan®) probes, labeled at one end with a fluorescent reporter dye (fluorophore) and at the other end with a quencher moiety. Probes labeled with different fluorophores are used to detect the various target sequences; up to 5 fluorophores may be used per reaction tube. The CTS 2.0 can detect a total of 27 DNA markers including 6 controls using 6 multiplex real-time PCR's in parallel. Control markers will assess the presence or absence of *Salmonella* and detect reaction failure. The absence or presence of each of the remaining 21 DNA markers generates a unique "genetic fingerprint", which is used to determine *Salmonella* type.

The CTS 2.0 assay uses dedicated software for result calling in a 2-step process. In the first step absence and presence of each of the 27 DNA markers is established. This generates a unique presence/absence pattern for all DNA markers, i.e., the "genetic fingerprint". In a second step this "genetic fingerprint" is matched to a database containing the genetic fingerprints of many *Salmonella* serovars. If a match is found the software will call a serovar result from the CTS 2.0 assays; if no match is found but *Salmonella* genus is confirmed a genovar result will be called.

Kit contents	Quantity
CTS 2.0 PCR Reagent Strips <i>PCR strips (with caps) having 12 tubes each, and containing 2 times 6 sets of target-specific primers and fluorescent probes in dry format</i>	3 pouches each having 8 PCR reagent strips with 12 tubes
CTS 2.0 PCR Master Mix <i>PCR Master Mix with DNA polymerase, dNTP's, reaction buffer and stabilizers</i>	2 tubes with 1 mL each
CTS 2.0 Chimney Caps <i>For 0.2 mL block only. Caps of the PCR strips need to be replaced with Chimney Caps when using 0.2 mL Block</i>	25 Chimney Caps
Sample Buffer	1 vial with 12 mL
Colony Samplers	90

3 EQUIPMENT AND MATERIALS REQUIRED BUT NOT PROVIDED

- ThermoFisher QuantStudio™ 5 Real-Time PCR System (96-well 0.1 mL or 0.2 mL Block)
- Vortex Mixer
- Heating Block for Eppendorf Tubes
- Centrifuge for Eppendorf Tubes
- Centrifuge for microtiter plates
- Holder for Eppendorf Tubes
- Eppendorf Tubes
- Holder for PCR Reagent Strips
- Pipettes & disposable hydrophobic filter tips for volumes of 20 – 200 µL
- Lab coat and powderless disposable gloves

4 WARNINGS AND PRECAUTIONS

- Caution - This procedure uses / detects pathogenic microorganisms and / or their metabolic products. Care should be taken to avoid ingestion or inhalation of potentially infectious aerosols, or contact with the skin. Laboratory personnel should follow appropriate laboratory safety precautions and have ready access to associated material safety data sheets (SDS, <http://www.ilpi.com/msds/>) While handling pathogens, personnel should use appropriate biosafety containment (<https://www.cdc.gov/labs/BMBL.html>) and be provided with appropriate personal protective equipment including clothing, protective gloves and appropriate eye protection
- The CTS 2.0 assay is for laboratory use only.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective pouches are open or broken upon arrival.
- Close protective pouches of reagents promptly with the zip seal after each use. Remove any excess air in the pouches prior to sealing.
- Do not remove desiccant from reagent pouches.
- Do not use reagents if desiccant is not present or is broken inside reagent pouches.
- Do not use reagents if the foil has been broken or damaged.
- Do not mix reagents from different pouches and/or kits and/or lots.
- Do not use expired reagents and/or materials.
- Do not expose the PCR Reagent Strips to excessive light to prevent degradation of fluorophores by photo-bleaching.
- Use powder-free gloves and avoid fingerprints and writing on PCR Reagent Strip caps.
- Good laboratory technique is essential to the proper performance of this assay. Due to the high analytical sensitivity of this test, extreme care should be taken to preserve the purity of all materials and reagents.
- To avoid contamination by amplicons, do not open the CTS 2.0 PCR Reagent Strips after use, but dispose immediately.
- Always handle specimens in accordance with safe laboratory procedures for Microbiological Laboratories.
- Wear protective clothing and disposable gloves while handling all reagents.
- Wash your hands thoroughly after performing the test.
- Do not smoke, drink, chew or eat in areas where specimens or kit reagents are being handled.
- Dispose of unused reagents and waste in accordance with local, state, provincial and/or federal regulations.

5 STORAGE AND STABILITY

CTS 2.0 Assay kit is stable at 2–8 °C through the stated expiration date. Do not use expired components. CTS 2.0 PCR Reagent Strips are provided in sealed pouches. To protect the product from humidity, immediately re-seal after opening.

PCR Reagents Strips are stable for up to 30 days at 2–8 °C after initial opening and re-sealing of the pouch.

6 INSTRUCTIONS FOR USE

Generate Work List

1. Log in to Check-Points CTS 2.0 Analysis portal with e-mail address and password.
2. Open the menu by clicking the  icon on the top left, when not visible.
3. Open the run sessions tab by clicking “Run sessions”.
4. Add a new run session with clicking the  icon on the top right. Global Identifier is automatically generated and cannot be changed. Custom run identifier may be entered under “Details”. Select the “Lot number” from the drop-down menu and select the correct “qPCR cyclers” to use (both mandatory). “Custom run identifier” and “Run session notes” are optional. NOTE: if QuantStudio™ 5 0.2 mL block is selected a pop-up will appear that the standard caps of the strips will need to be replaced with the supplied Chimney Caps. The standard caps can be used on the 0.1mL block To add samples, click the pipette icon  and enter sample information. “Samples names” are mandatory, other information like ‘Sample location’, Production batch”, Sample origin”, “Sampling device” and “Sampling date” may also be entered.
5. Click “Save + New” to go to the next sample and continue until the final sample and then click “Save”.
6. Alternatively, click the import samples icon  to upload a CSV or XLSX file with sample information. File format may be downloaded using the export button .
7. Click “Download Cycler run file”. This generates a plate setup for the real-time PCR run in the “PCR System Operation” step. If sample info is changed afterwards, always click “Download Cycler run file” again to have the most recent .txt file for the cycler. Running the test with the outdated .txt plate setup will result in an error during analysis
8. Next, download and print a work list for use in the laboratory by clicking the  icon on the top right.
9. Proceed to Sample preparation.

Sample preparation

NOTE: One (1) PCR Reagent Strip is used for two (2) samples. In case of uneven sample numbers supplement with samples from a control (reference) strain.

1. Turn on heating block at $98^{\circ}\text{C} \pm 3^{\circ}$.
2. Dispense **200 μL** Sample Buffer into a 1.5 mL tube. Use a separate tube for each sample and write sample ID on the tube.
3. Pierce through a single colony into the agar using a colony sampler. Briefly touch the bottom of the plate and take the sampler out again. Always keep the colony sampler in a vertical position as shown in Figure 1.




Figure 1: Correct colony sampling

4. Place the colony sampler in 1.5 mL tube with **200 μL** of Sample Buffer and roll the sampler between thumb and index finger at least 5 times back and forth while remaining in the buffer.
5. Remove and discard the colony sampler.
6. Close the tube and vortex briefly.
7. Place tubes in pre-heated heating block and incubate for **10 min at $98^{\circ}\text{C} \pm 3^{\circ}$** .
8. Vortex tubes and place them on laboratory bench to cool down to room temperature for at least 5 minutes. If the samples are not used directly spin down briefly and store tubes at -20°C .
9. Proceed to “Real-time PCR Set-Up”.

Real-time PCR Set-Up

NOTE: Samples from -20°C should be placed at room temperature for 15 minutes, vortexed and spun down briefly.

1. Dispense **35 µL** PCR Master Mix into a 1.5 mL tube. Use a separate tube for each sample and write sample ID on the tube.
2. Briefly spin down the samples from the previous step (Sample preparation), and add **140 µL** sample to the PCR Master Mix. Mix by pipetting up and down gently 3 times.
3. Take the required number of PCR Reagent Strips from their protective pouches. Remove excess air, and close pouches with the zip seal. Label the end of each PCR Reagent Strip with the appropriate specimen identification.
4. Remove strip caps carefully and keep separate.




Figure 2 Pipetting scheme of Sample/Master mix specimens.

5. Transfer **25 µL** of the first Sample/Master Mix specimen to tube 1 of PCR Reagent Strip (Figure 2). Repeat this for tubes 2 – 6.
6. Transfer **25 µL** of the second Sample/Master Mix specimen to tube 7 its PCR Reagent Strip. Repeat this for tubes 8 – 12 (Figure 2).
7. Proceed until all Sample/Master Mix specimens have been transferred to their PCR Reagent Strips.
8. Carefully place the strip caps back on the strips. Make sure strip caps are sealed airtight to prevent evaporation during PCR, **and spin down briefly.**
IMPORTANT: if a 0.2 mL block is used, replace the standard caps with the supplied chimney caps to close the strips.
9. Power on the PCR System and computer.
10. Open the drawer of the real-time PCR cycler.
11. Place PCR Reagent Strips in the PCR machine using the worklist created in “Generate work list” section. If needed place “dummy” PCR reagent strips according to the worklist to distribute strips evenly over the PCR heating block.
12. Close the drawer of the real-time PCR cycler.
13. Proceed to “PCR System Operation”.

PCR System Operation (DA software version 1.x)

NOTE: It is recommended to use an instrument connected to a computer or laptop. In absence of a direct connection please follow the alternative steps described in the IFU.

1. Start the QuantStudio™ Design & Analysis Software (1.x).
2. Select **"File"** go to **"Open..."** and select the template file Q5 Xml CTS 2 RUN FILE.edt (X= 0_1 mL or 0_2 mL).
3. Select **"File"** go to, **"Import Plate Setup..."** click **"Browse"** and select the .txt file generated by portal and select **"Apply"**. If the import of the file was successful, a pop-up with Imported successfully will be shown.
4. Select OK.
5. Continue with step 6 if the computer is attached to the cycler, otherwise continue with step 5.1
 - 5.1. Select **"File"** and select **"Save As..."** and save the file to a USB stick, changing the name to the same name as the downloaded .txt cycler run file from the portal (copy-paste name from .txt file without the .txt).
Proceed with step 6.1
6. Select **"Start Run"** in the Run tab and select the correct machine to initiate the run. Save the run without changing the name of the file in the pop-up. Continue with step 7
 - 6.1. When working without computer attached to cycler: Plug the USB stick into the PCR system. On the touch screen, select **"Load Experiment"**, **"USB"**, select the file created under 5.1, and **"Start Run"**
7. At the end of the run, save (or transfer) the data file (.eds extension) to a dedicated network location or USB device
8. Make sure the file from step 7 is (still) open in the software, go to the Export tab, and select under File type "RDML" (make sure "Open exported files when complete" is deselected)
Select the dedicated network location or USB device as the location you want the file to be saved to using the Browse button, and click **"Export"** (do not change the name)

PCR System Operation (Design and Analysis software DA2.x and DA3.x)

1. Start the QuantStudio™ Design & Analysis Software (2.x or 3.x)
2. On the software dashboard, select "Open file" located near the top left side (DA2.x) or right side (DA3.x) of the window.
3. Select the template file Q5 Xml CTS 2 RUN FILE.edt, depending on the instrument's block type) (X= 0_1 mL or 0_2 mL).
4. Proceed to the plate setup tab. Select the three dots next to the gear wheel, located above the plate map. Select "Import plate Setup" and chose the downloaded .txt cycler run file from the portal. Click 'OK' when prompted to overwrite the template plate setup.
5. On the right-hand side of the window, select "Actions -> save as" and save the template file (.edt extension) using the identical run file name generated by the CTS2.0 portal (copy-paste name from .txt file without the .txt extension).

IMPORTANT: Incorrect file names will lead to errors during analysis of exported data through the Check-Points cloud software. If upload issues occur, please contact CTS for support.

Continue with step 8 if the computer is attached to the QuantStudio™ 5 instrument, otherwise continue with step 6+7

6. Transfer the saved template file (.edt extension) from step 5 to a USB stick. Plug USB stick in the PCR system. On the touch screen of the QS5, select **"Load Experiment"**, **"USB"**, select the file created in step 5 and **"Start Run"**
7. At the end of the run save and transfer the data file (.eds extension) from the instrument to a dedicated network location or USB device. And continue with step 9.
8. Select **"Start Run"** in the Run tab and select the correct machine to initiate the run. Save the run without changing the name of the file in the pop-up.

9. Ensure the file from step 8 is (still) open in the DA software. Otherwise open the .eds file transferred from step 7. Perform Analyse if necessary, then select “Actions -> export to RDML”. Select the dedicated network location or USB device as the location you want the file to be saved to using the Browse button and click “**Export**”, do not change the name (check that the prompted file name is the same as the downloaded .txt cyclor run file from the portal).

Data Processing and Results

1. Log in to the Check-Points CTS 2.0 Analysis portal.
2. Open the Run session generated in the “Generate work list” section.
3. Click the “**Upload Cyclor Result**” button and upload the exported RDML file containing the raw PCR data.
4. The data processing runs automatically in the background and it is not necessary to keep the window open. Processing will take a few minutes and results will pop-up automatically. If results do not appear, refresh the browser once (F5) and check the “**Analysis**” tab.
5. Reports can be printed by clicking the “**Reporting**” tab. Use the tick boxes to select results for printing and use  to download individual reports or  to download a summary report.
6. Table 1 below shows the results and interpretations that may be obtained.

Table 1 CTS Assay Result Interpretation

Assay Result Reported	Interpretation of Result	Next Step
Salmonella Serovar	Salmonella positive and Serovar result	Final result
Salmonella Serovar 1 or Serovar 2	Salmonella positive result with two possible Serovar results	Final result
Salmonella positive Genovar	Salmonella positive result with Genovar code	Repeat sample once with fresh isolate
Salmonella positive, Inconclusive typing result	Salmonella positive, but no typing result	Repeat sample once with fresh isolate
No Salmonella	Gram-negative bacteria detected but not Salmonella	Final result
No DNA Detected	No gram-negative bacteria detected	Repeat sample once with fresh isolate
Reaction not OK	Reaction failed – inhibitory sample or reagent failure; no Internal Control amplification	Repeat sample once with fresh isolate

Repeat procedure

In the following circumstances it is advised to repeat the sample once with a fresh isolate:

1. Salmonella positive with Genovar code
2. Salmonella positive with inconclusive typing result
3. No DNA detected
4. Reaction not OK

A Salmonella positive with Genovar code or Inconclusive typing result may be obtained in case one (1) of the 21 typing markers has given an erroneous result generating a Genovar code, for which there is no association with a known Serovar in the CTS 2.0 database. This is a rare event for individual markers, but significant for the sum of the 21 typing markers. Repeating the sample once will generate a Serovar result if this specific issue has occurred and should be noted as the final result. The initial Genovar result or Inconclusive typing result may also be the proper result and repeating the assay will generate the same result in these cases.

In case of "Reaction not OK" or "No DNA detected" it is also advised to repeat the sample once with a fresh isolate. These results are often caused by sub-optimal reaction conditions and/or impurities in the sample. Repeating the reaction with a fresh isolate will generate a proper result the second time in many cases.

It may be that the Check&Trace Salmonella 2.0 Assay does not generate a conclusive confirmation and/or typing result, even after repeating the test. In such cases it is advised to use an alternative method to obtain a final test result. (1,2).

Requesting support

Support about the obtained results can be requested by marking the support button in the Run session.

1. Select the run session for which you like to request support
2. Select the "person" icon.
3. Click confirm on the pop-up




Figure 3 Image showing the "person" support request icon (yellow circle) and the pop-up after clicking.

Figure 4 Message center after clicking the support button

4. **Add** a descriptive **message** so support knows where to help. Click save to finalize
5. Adding a message is **mandatory**, otherwise a support request won't be sent by the system.

7 QUALITY CONTROL

Quality control procedures monitor the performance of the assay. Laboratories must establish the number, type and frequency of testing control materials according to guidelines or requirements of local, provincial, state and/or federal regulations or accreditation organizations in order to monitor the entire analytical process.

External Positive and Negative Controls are not used by the CTS 2.0 Assay for the purpose of test result interpretation. The CTS 2.0 Assay also has a total of 6 internal controls to monitor for potential reaction failure. However, it is advised to test control samples daily until adequate process validation is achieved. Reduced frequency of control testing should be in accordance with applicable regulations.

External Controls should yield the expected results: Salmonella serovar for External Positive Control, and "No Salmonella" for External Negative Control. Any other result for the external controls is indicative for system failure.

Table 2 below specifies a number of reference strains that may be used for as external positive or negative controls. Other controls strains of which the test result is well documented may also be used.

Table 2 Commercially Available Strains for External Positive and negative Control

External Control Strain	Test Result
<i>Salmonella enterica subsp. enterica serovar</i> Enteritidis (NCTC 12694)	Salmonella Enteritidis
<i>Salmonella enterica subsp. enterica serovar</i> Typhimurium (NCTC 10413)	Salmonella Typhimurium
<i>Salmonella enterica subsp. enterica serovar</i> Typhimurium (Monophasic); I 4,5,12:i:- (NCTC 13952)	Salmonella Typhimurium Monophasic
<i>Escherichia coli</i> (ATCC 25922)	No Salmonella


8 LIMITATIONS OF THE PROCEDURE

- This product can only be used on the ThermoFisher QuantStudio™ 5 Real-Time PCR System (96-well 0.1 mL or 0.2 mL Block)
- The test is not intended for IVD or human/animal health testing.
- Erroneous results may occur from improper sample collection, handling, storage, technical error, sample mix-up.
- If the CTS 2.0 result is "*Salmonella* positive with genovar code", "*Salmonella* positive inconclusive typing result", "No DNA detected" or "Reaction not OK", the test should be repeated with a fresh isolate.
- The following serovars can be detected:
Abaetetuba, Aberdeen, Abortusovis, Adelaide, Agama, Agona, Alachua, Albany, Altona, Amsterdam, Anatum, Banana, Bareilly, Berta, Blockey, Bovismorbificans, Bracknell, Braenderup, Brandenburg, Bredeney, Cerro, Chailey, Chester Choleraesuis, Corvallis, Cubana, Derby, Dublin, Eastbourne, Emek, Enteritidis, Fresno, Gallinarum Gallinarum, Gallinarum Pullorum, Gaminara, Give, Gloucester, Goldcoast, Hadar, Haifa, Hartford, Havana, Heidelberg, Hvitittingfoss, Idikan, Indiana, Infantis, Isangi, Javiana, Jerusalem, Johannesburg, Kapemba, Kedougou, Kentucky, Kiambu, Kingston, Kottbus, Lexington, Litchfield, Liverpool, Livingstone, Llandoff, London, Manhattan, Mbandaka, Meleagridis, Minnesota, Molade, Montevideo, Muenchen, Muenster, Newport, Nottingham, Offa, Ohio, Oranienburg, Ordonez, Orion, Ouakam, Panama, Paratyphi B biovar. Java, Poona, Putten, Reading, Rissen, Saintpaul, Salmonella Bongori, Sandiego, Schwarzengrund, Senftenberg, Soerenga, Stanley, Taksony, Tennessee, Thompson, Typhimurium, Uganda, Urbana, Virchow, Weltevreden, Worthington, Yoruba, monophasic variant of Salmonella Typhimurium (1,4,[5],12:i:-), 4,[5],12:d:-, 47:z4,z23:-, 4,12:i:- (not Typhimurium)
- *Salmonella* serovars other than the above 106 types will yield a genovar code.
- In some cases, CTS 2.0 may report "*Salmonella* positive serovar A or serovar B".
- Mutations or polymorphisms in primer or probe binding regions may affect detection of certain sequence variants (Sequence Types) of the serovars included in the CTS 2.0 assay. In such cases CTS 2.0 will report a genovar code instead of the serovar result.
- CTS 2.0 is able to discriminate between *Salmonella enterica* and *Salmonella bongori*, but not between the various *enterica* subspecies *arizonae*, *diarizonae*, *enterica*, *houtenae*, *indica* and *salamae*. In all these cases the CTS 2.0 result will be "*Salmonella* positive genovar code".
- CTS 2.0 is a qualitative test and does not provide quantitative values nor indicate the quantity of organisms present.
- The performance of the CTS 2.0 has been evaluated on Nutrient Agar and XLD agar. Other media will need to be evaluated separately for any laboratory to verify adequate performance on such media.

9 REFERENCES

1. ISO 6579-1. Microbiology of the food chain – Horizontal method for detection, enumeration and serotyping of *Salmonella* - Part 1: Detection of *Salmonella* spp. [ISO 6579-1:2017].
2. ISO 6579-3. Microbiology of the food chain – Horizontal method for detection, enumeration and serotyping of *Salmonella* - Part 3: Guidelines for serotyping of *Salmonella* spp. [ISO 6579-3:2014].
3. Antigenic formulae of the *Salmonella* serovars. Grimont, P.A.D. and Weill, F.X. Vol.9, 2007.

10 SYMBOLS GLOSSARY

Symbol	Meaning
	Manufacturer
	Contains sufficient for <n> tests
	Temperature limit
	Consult instructions for use For electronic Instructions for use, the URL accompanies the symbol.
	Keep dry
	Keep away from light

