

OvaCyte™

Welcome to the future of parasitology

OvaCyte is a compact and novel point-of-care automated system that uses digital analysis and Artificial Intelligence (AI) to count and identify eggs of intestinal parasites and trematodes in a wide range of host species.



**Ultrafast, easy,
smart & reliable**

Sample preparation takes less than 3 minutes of hands-on time. A fully automated process then delivers fast, reliable results. No technical training required.



**AI-powered
Software**

OvaCyte's AI-powered software eliminates manual microscopy, leveraging a library of over 1 million images and 150,000 tagged parasites.



**Accurate results
in minutes**

Results are available in real time on the device and immediately after the scan in the OvaCyte Web App together with all scan images and reports.



MORE INSIGHTS WITH OVINE SPECIATION

Species- and group-level differentiation, including the % of highly pathogenic coccidia and *Haemonchus contortus*.

TEST TWO SAMPLES AT ONCE

Save time and 66% plastic with OvaCyte DUAL Cassettes.



Coverage of key parasites

Standard FEC test

(OvaCyte EQ & LA FEC PLUS Test Kit)

Equine	Bovine/Ovine
Strongyles	Strongyles
<i>Parascaris</i> spp.	Coccidia (<i>Eimeria</i> spp.)
<i>Strongyloides westeri</i>	<i>Moniezia</i> spp.
<i>Anoplocephala</i> spp.	<i>Nematodirus</i> spp.
<i>Oxyuris equi</i>	<i>Strongyloides papillosus</i> spp.

OvaCyte Speciation Ovine

(Optional software feature, activated via credits)

Ovine
Identification of <i>Haemonchus contortus</i>
Percentage of suspected highly pathogenic coccidia (<i>Eimeria ovinoidalis</i> , <i>E. crandallis</i> , and <i>E. weybridgeensis</i>)
Identification of <i>Nematodirus battus</i>

Fluke test

(OvaCyte LA Fluke PLUS Test Kit)

Bovine/Ovine
<i>Fasciola hepatica</i>
<i>Paramphistomum</i> spp.
NEW! Equine test (OvaCyte EQ Ano Test Kit)
<i>Anoplocephala</i> spp.
<i>Oxyuris equi</i>

OvaCyte is as easy as A, B, C.



Collect a fecal sample and prepare it using the appropriate solutions, following the instructions specific to the test kit being used. Once the sample is prepared, fill the cassette according to the kit instructions, then place it on the device to start the analysis. After a short analysis time, the results are automatically generated and can be viewed on the OvaCyte device as well as via the web app on a computer or smartphone.

Ordering information

Product*	Short description/contents	Cat. No.
OvaCyte Fecal Analyzer	Compact device with integrated touch screen, including power cords, power adapter, locking nut and Ethernet cable	TN-FP1001-02A For U.S.: FP1001-USA
Universal EQ & LA Starter Pack with Fluke	Everything you need to kickstart your journey with OvaCyte EQ & LA and Fluke	FP25001
OvaCyte EQ & LA FEC PLUS Test Kit	40 x EQ & LA Cassettes and 2 x 1 liter of EQ & LA Flotation Fluid	FP27001 (40)
OvaCyte DUAL EQ & LA Test Kit	20 x Dual EQ & LA Cassettes and 2 x 1 liter of EQ & LA Flotation Fluid	FP42001 (40)
OvaCyte Speciation Ovine	Precision AI to estimate the percentage of <i>Haemonchus contortus</i> among strongyles, suspected highly pathogenic coccidia, and <i>Nematodirus battus</i> among <i>Nematodirus</i> spp. in sheep fecal samples. Activated via credits, compatible with the PLUS and DUAL test kits.	SP0001 (10) SP0002 (50) SP0003 (100)
NEW! OvaCyte EQ Ano Test Kit	40 x EQ Ano Cassettes and 750 ml of EQ Ano Flotation Fluid	FP52001 (40)
OvaCyte LA Fluke PLUS Test Kit	40 x Fluke Cassettes and 1 liter of Fluke Flotation Fluid	FP16001 (40)

*Products are not intended for the diagnosis, prevention or treatment of a disease. Two syringes are provided with each kit.



Scan QR code to see
OvaCyte in action



Scan QR code to learn more about
OvaCyte and book a free demo

Get in touch with our specialists: www.indical.com/ovacyte | Request a quote or order online: shop.indical.com

Email (Americas): us_support@indical.com | us_orders@indical.com | Call (Americas): **+1 (407) 766 5128**

Email: support@indical.com | orders@indical.com | Call: **+49 341 124 54 0**

Connect with us on LinkedIn: www.linkedin.com/company/indical

Products are not intended for the diagnosis, prevention or treatment of a disease. For up-to-date licensing information and product-specific disclaimers, see the respective handbook or user manual. Regulatory requirements vary by country, products may not be available in your geographic area. Trademarks: INDICAL® (INDICAL BIOSCIENCE GmbH), OvaCyte™ (Telenostic Ltd.). Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are not to be considered unprotected by law. Product images may differ from the actual product. INFLY-OvaCyte-EN-261006 © 2026 INDICAL BIOSCIENCE

INDICAL
BIOSCIENCE

INDICAL BIOSCIENCE GmbH
Deutscher Platz 5b
04103 Leipzig, Germany