

NEW! IndiMix JOE

Build your assays faster and without compromise

- **Novel master mix:** Built-in, validated internal control for amplification OR extraction
- **Maximum flexibility:** One protocol for DNA and RNA applications
- **Compatible:** Works with common PCR thermocyclers
- **High-quality performance:** Inhibitor tolerant, reproducible results, even in challenging applications



One master mix. One protocol. For results you can trust.
Accelerate your assay development with our novel IndiMix JOE.

A master mix compatible with your laboratory

Figure 1



IndiMix JOE is compatible with common PCR thermocyclers and can be used to identify even the smallest number of copies per well. Get great results with both fast and standard mode cyclers.

Figure 1. A dilution series of *in vitro* RNA shows reliable performance with common PCR cyclers, down to 10 copies/well.

Extensive coverage of sample types and targets

Figure 2

Sample	DNA C _t	RNA C _t
Lung A	Negative	18.71
Lung B	7.49	Negative
Lymph node A	16.74	Negative
Serum A	20.78	Negative
Serum B	20.72	Negative
Nasal swab A	Negative	32.28
Nasal swab B	Negative	33.09

IndiMix JOE was developed and tested to work with a wide variety of animal sample matrices and applications.

INDICAL's novel IndiMix JOE master mix enables real-time amplification of single or multiple nucleic acid targets and contains both reverse transcriptase and polymerase to amplify DNA and RNA targets.

Figure 2. DNA and RNA targets from different sample types tested.

Request a quote or order online: shop.indical.com

Email: support@indical.com | orders@indical.com | Phone: +49 341 124 54 0 | Fax: +49 341 124 54 60

Email (Americas): us_support@indical.com | us_orders@indical.com | Phone (Americas): +1 (407) 766 5128

Connect with us on LinkedIn: www.linkedin.com/company/indical

For up-to-date licensing information and product-specific disclaimers, see the respective handbook or user manual. Regulatory requirements vary by country, products may not be available in your geographic area. For molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease. Trademarks: INDICAL®, IndiMix™ (INDICAL BIOSCIENCE GmbH). Bio-Rad® (Bio-Rad Laboratories, Inc.). JOE™ (Life Technologies Corporation). Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are not to be considered unprotected by law. Product images may differ from the actual product. INFLY-IndiMix-EN-202201 © 2022 INDICAL BIOSCIENCE

INDICAL
BIOSCIENCE

INDICAL BIOSCIENCE GmbH
Deutscher Platz 5b
O4103 Leipzig, Germany

INDICAL Inc.
1317 Edgewater Dr #3722
Orlando, FL 32804, USA

