

IndiMix JOE & IndiMix TAMRA

Build your assays faster and without compromise


- **Novel master mix:** Built-in primers and probes to detect extraction or amplification internal control
- **Maximum flexibility:** One protocol for DNA and RNA applications – Choose from JOE or TAMRA as internal control label
- **Efficiency:** PCR results in less than 1 hour
- **Robust performance:** Inhibitor tolerant and reproducible results, including challenging applications



Reliable results from one master mix and one protocol
Accelerate assay development with IndiMix

A master mix compatible with any laboratory

Data shown for IndiMix JOE



IndiMix master mixes are compatible with common real-time PCR thermocyclers and can be used to identify a minimal number of copies per reaction. Obtain accurate results with both fast and standard mode cyclers.

Figure 1. A dilution series of *in vitro* RNA shows reliable performance of IndiMix JOE with common PCR cyclers, down to 10 copies/well.

Amplification for various sample types and targets

Data shown for IndiMix JOE

Sample	DNA C _T	RNA C _T
Lung A	Negative	18.71
Lung B	7.49	Negative
Lymph node A	16.74	Negative
Serum A	20.78	Negative
Serum B	20.72	Negative
Nasal swab A	Negative	32.28
Nasal swab B	Negative	33.09

IndiMix products were developed and tested with a multitude of biological sample matrices. IndiMix master mixes enable real-time amplification of single or multiple nucleic acid targets and contain both polymerase and reverse transcriptase to amplify DNA and RNA targets.

Figure 2. Representation of samples tested for DNA and RNA with IndiMix JOE.



Scan QR code to learn more about IndiMix TAMRA



Scan QR code to learn more about IndiMix JOE

Request a quote or order online: shop.indical.com

Email (Americas): us_support@indical.com | us_orders@indical.com | Call (Americas): +1 (407) 766 5128

Email: support@indical.com | orders@indical.com | Call: +49 341 124 54 0

Connect with us on LinkedIn: www.linkedin.com/company/indical

For up-to-date licensing information and product-specific disclaimers, see the respective handbook or user manual. Regulatory requirements vary by country, products may not be available in your geographic area. For molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease. Trademarks: INDICAL®, IndiMix™ (INDICAL BIOSCIENCE GmbH). Bio-Rad® (Bio-Rad Laboratories, Inc.). JOE™, TAMRA™ (Thermo Fisher Scientific). Registered names, trademarks, etc., used in this document, even when not specifically marked as such, are not to be considered unprotected by law. Product images may differ from the actual product. INFLY-IndiMix-JOE-TAMRA-EN-231218 © 2023 INDICAL BIOSCIENCE

INDICAL
BIOSCIENCE

INDICAL BIOSCIENCE GmbH
Deutscher Platz 5b
O4103 Leipzig, Germany

INDICAL Inc.
1317 Edgewater Dr #3722
Orlando, FL 32804, USA

