

**Trade name:** Buffer MVL

**Current version :** 1.0.1, issued: 06.08.2025

**Replaced version:** 1.0.0, issued: 02.02.2021

**Region:** GB

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1 Product identifier

**Trade name**

**Buffer MVL**

**Form**

The product is included in the following product sets:

**SP69805 IndiSpin Mastitis Kit (192)**

**SP947757 IndiMag Mastitis Kit (384)**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses of the substance or mixture**

laboratory chemical

**Uses advised against**

No data available.

### 1.3 Details of the supplier of the safety data sheet

**Address**

INDICAL BIOSCIENCE GmbH

Deutscher Platz 5b

04103 Leipzig

Germany

Telephone no. +49 341 12454 0

Fax no. +49 341 12454 60

e-mail [compliance@indical.com](mailto:compliance@indical.com)

**Advice on Safety Data Sheet**

[sdb\\_info@umco.de](mailto:sdb_info@umco.de)

### 1.4 Emergency telephone number

For medical advice (in German and English):

+49 (0)551 192 40 (Giftinformationszentrum Nord)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

**Classification in accordance with Regulation (EC) No 1272/2008 (CLP)**

Acute Tox. 4; H302

Acute Tox. 4; H332

Aquatic Chronic 3; H412

Eye Dam. 1; H318

Skin Corr. 1C; H314

**Classification information**

This product is assessed and classified using the methods and criteria below referred to in Article 9 of Regulation (EC) n° 1272/2008:

Physical hazards: determined through assessment data based on the methods or standards referred to in part 2 of Annex I to CLP

Health hazards and environmental hazards: determined through toxicological and ecotoxicological assessment data based on the methods or standards referred to in Part 3, 4 and 5 of Annex I to CLP.

### 2.2 Label elements

**Labelling according to Regulation (EC) No 1272/2008 (CLP Regulation)**

**Hazard pictograms**



Trade name: Buffer MVL

Current version : 1.0.1, issued: 06.08.2025

Replaced version: 1.0.0, issued: 02.02.2021

Region: GB

GHS05

GHS07

**Signal word**

Danger

**Hazardous component(s) to be indicated on label:**

GUANIDINIUMTHIOCYANATE

**Hazard statement(s)**

H302+H332

Harmful if swallowed or if inhaled

H314

Causes severe skin burns and eye damage.

H412

Harmful to aquatic life with long lasting effects.

**Hazard statements (EU)**

EUH032

Contact with acids liberates very toxic gas.

**Precautionary statement(s)**

P260

Do not breathe mist/vapours/spray.

P280

Wear protective gloves/protective clothing/eye protection/face protection.

P303+P361+P353

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].

P305+P351+P338

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310

Immediately call a POISON CENTER/doctor.

**2.3 Other hazards**

PBT assessment

The product is not considered to be a PBT.

vPvB assessment

The product is not considered to be a vPvB.

**SECTION 3: Composition/information on ingredients****3.1 Substances**

Not applicable. The product is not a substance.

**3.2 Mixtures****Hazardous ingredients**

No	Substance name		Additional information	
	CAS / EC / Index / REACH no	Classification (EC) 1272/2008 (CLP)	Concentration	%
1	<b>GUANIDINIUMTHIOCYANATE</b>			
	593-84-0	Acute Tox. 4*; H302	>= 25.00 - < 50.00	wt%
	209-812-1	Acute Tox. 4*; H312		
	615-004-00-3	Acute Tox. 4*; H332		
	-	Skin Corr. 1C; H314		
	-	Eye Dam. 1; H318		
	-	Aquatic Chronic 3; H412		
	-	EUH032		

Full text of H- and EUH-phrases, if not already mentioned in section 2.2: see section 16.

(\*, \*\*, \*\*\*, \*\*\*\*) Detailed explanation pls. refer to CLP regulation No. 1272/2008, annex VI, 1.2

**SECTION 4: First aid measures****4.1 Description of first aid measures****General information**

Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing. Seek medical advice immediately. Show MSDS to the attending physician.

**After inhalation**

Remove affected persons from dangerous area by observing suitable respiratory protection measures. Ensure supply of fresh air. Do not use mouth-to-mouth or mouth-to-nose resuscitation. In case of persisting adverse effects consult a physician.

**Trade name:** Buffer MVL**Current version :** 1.0.1, issued: 06.08.2025**Replaced version:** 1.0.0, issued: 02.02.2021**Region:** GB**After skin contact**

Wash off immediately with soap and water. Seek medical attention.

**After eye contact**

Remove contact lenses. Rinse eye thoroughly under running water keeping eyelids wide open and protecting the unaffected eye (at least 10 to 15 minutes). Get immediate ophthalmic treatment.

**After ingestion**

Rinse the mouth thoroughly with water. Do not induce vomiting. Never give anything by mouth to an unconscious person.

**4.2 Most important symptoms and effects, both acute and delayed**

No data available.

**4.3 Indication of any immediate medical attention and special treatment needed**

No data available.

**SECTION 5: Firefighting measures****5.1 Extinguishing media****Suitable extinguishing media**

Water spray jet; Alcohol-resistant foam; Dry chemical extinguisher; Carbon dioxide

**Unsuitable extinguishing media**

High power water jet

**5.2 Special hazards arising from the substance or mixture**

In the event of fire, the following can be released: Carbon oxides (COx); Sulphur oxides (SxOy); Nitrogen oxides (NOx)

**5.3 Advice for firefighters**

Use self-contained breathing apparatus. Wear protective clothing. Run-off water from fire fighting must not be discharged into drains or enter surface water.

**SECTION 6: Accidental release measures****6.1 Personal precautions, protective equipment and emergency procedures****For non-emergency personnel**

Refer to protective measures listed in sections 7 and 8. Use personal protective clothing.

**For emergency responders**

Personal protective equipment (PPE) - see section 8.

**6.2 Environmental precautions**

Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil. In case of entry into waterways, soil or drains, inform the responsible authorities.

**6.3 Methods and material for containment and cleaning up**

Contain and collect spillage with non-combustible absorbent materials, e.g. sand, earth, vermiculite, diatomaceous earth and place in container for disposal according to local regulations (see section 13).

**6.4 Reference to other sections**

Information regarding safe handling, see section 7. Information regarding personal protective measures, see section 8. Information regarding waste disposal, see section 13.

**SECTION 7: Handling and storage****7.1 Precautions for safe handling****Advice on safe handling**

Risks inherent to handling the product must be minimised by applying the appropriate protective and preventive measures. Working processes should - so far as possible, according to the state of the art - be designed to rule out bodily contact or the release of hazardous substances.

**General protective and hygiene measures**

**Trade name:** Buffer MVL

**Current version :** 1.0.1, issued: 06.08.2025

**Replaced version:** 1.0.0, issued: 02.02.2021

**Region:** GB

Do not eat, drink or smoke during work time. Keep away from foodstuffs and beverages. Do not inhale vapours. Avoid contact with eyes and skin. Wash hands before breaks and after work. Remove contaminated clothing and shoes and launder thoroughly before reusing. Have emergency shower available. Provide eye wash fountain in work area.

**Advice on protection against fire and explosion**

Keep away from ignition sources and provide for good ventilation.

**7.2 Conditions for safe storage, including any incompatibilities**

**Technical measures and storage conditions**

Keep container tightly closed and dry in a cool, well-ventilated place.

**Requirements for storage rooms and vessels**

Containers which are opened must be carefully closed and kept upright to prevent leakage. Always keep in containers of same material as the original.

**Incompatible products**

Do not store together with: Acids

**7.3 Specific end use(s)**

No data available.

**SECTION 8: Exposure controls/personal protection**

**8.1 Control parameters**

No parameters available for monitoring.

**8.2 Exposure controls**

**Appropriate engineering controls**

Provide adequate ventilation. Where reasonably practicable this should be achieved by the use of local exhaust ventilation and good general extraction. If these are not sufficient to maintain concentrations of particulates and solvent vapour below the OEL (=Occupational Exposure Limit), suitable respiratory protection must be worn.

**Personal protective equipment**

**Respiratory protection**

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of aerosol and mist formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

**Eye / face protection**

Tightly fitting safety glasses (EN 166).

**Hand protection**

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product. Before use, the protective gloves should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

**Other**

Chemical-resistant work clothes. Protective shoes.

**Environmental exposure controls**

No data available.

**SECTION 9: Physical and chemical properties**

**9.1 Information on basic physical and chemical properties**

<b>State of aggregation</b>	liquid
<b>Form</b>	liquid
<b>Colour</b>	

Trade name: Buffer MVL

Current version : 1.0.1, issued: 06.08.2025

Replaced version: 1.0.0, issued: 02.02.2021

Region: GB

No data available

**Odour**

characteristic

**pH value**

Value 8

Source supplier

**Boiling point / boiling range**

No data available

**Melting point/freezing point**

No data available

**Decomposition temperature**

No data available

**Flash point**

No data available

**Ignition temperature**

No data available

**Auto-ignition temperature**

not determined

Source supplier

**Flammability**

No data available

**Lower explosion limit**

No data available

**Upper explosion limit**

No data available

**Vapour pressure**

No data available

**Relative vapour density**

No data available

**Relative density**

No data available

**Density**Value 1.14 g/cm<sup>3</sup>

Source supplier

**Solubility**

No data available

**Partition coefficient n-octanol/water (log value)**

No data available

**Kinematic viscosity**

No data available

**Particle characteristics**

No data available

**9.2 Other information****Other information**

No data available.

**SECTION 10: Stability and reactivity**

Trade name: Buffer MVL

Current version : 1.0.1, issued: 06.08.2025

Replaced version: 1.0.0, issued: 02.02.2021

Region: GB

**10.1 Reactivity**

Dangerous reactions are not expected if the product is handled according to its intended use.

**10.2 Chemical stability**

Stable under recommended storage and handling conditions (See section 7).

**10.3 Possibility of hazardous reactions**

Dangerous reactions are not to be expected when handling product according to its intended use.

**10.4 Conditions to avoid**

None, if handled according to intended use.

**10.5 Incompatible materials**

strong acids; Sodium hypochlorite; Oxidizing agents; Alkalis

**10.6 Hazardous decomposition products**

None, if handled according to intended use.

**SECTION 11: Toxicological information****11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**

<b>Acute oral toxicity (result of the ATE calculation for the mixture)</b>	
<b>Product Name</b>	
<b>Buffer MVL</b>	
ATE (Mixture)	1020.41 mg/kg
Method	Calculation method according Regulation (EC) No 1272/2008, (CLP), annex I, part 3, section 3.1.3.6.

<b>Acute oral toxicity</b>	
No data available	

<b>Acute dermal toxicity (result of the ATE calculation for the mixture)</b>	
<b>Product Name</b>	
<b>Buffer MVL</b>	
Comments	The result of the applied calculation method according to the European Regulation (EC) 1272/2008 (CLP), Paragraph 3.1.3.6, Part 3 of Annex I is outside the values that imply a classification / labelling of this mixture according to table 3.1.1 defining the respective categories (ATE dermal > 2000 mg/kg).

<b>Acute dermal toxicity</b>	
No data available	

<b>Acute inhalational toxicity (result of the ATE calculation for the mixture)</b>	
<b>Product Name</b>	
<b>Buffer MVL</b>	
ATE (Mixture)	3.0612 mg/l
Route of exposure / physical form	Dust/mist
Method	Calculation method according Regulation (EC) No 1272/2008, (CLP), annex I, part 3, section 3.1.3.6.

<b>Acute inhalational toxicity</b>	
No data available	

<b>Skin corrosion/irritation</b>	
<b>Product Name</b>	
<b>Buffer MVL</b>	
Source	supplier
Evaluation	strongly corrosive

<b>Serious eye damage/irritation</b>	
<b>Product Name</b>	

**Trade name:** Buffer MVL

**Current version :** 1.0.1, issued: 06.08.2025

**Replaced version:** 1.0.0, issued: 02.02.2021

**Region:** GB

<b>Buffer MVL</b>	
Source	supplier
Evaluation	Irreversible effects on the eye
<b>Respiratory or skin sensitisation</b>	
No data available	
<b>Germ cell mutagenicity</b>	
No data available	
<b>Reproduction toxicity</b>	
No data available	
<b>Carcinogenicity</b>	
No data available	
<b>STOT - single exposure</b>	
No data available	
<b>STOT - repeated exposure</b>	
No data available	
<b>Aspiration hazard</b>	
No data available	
<b>Endocrine disrupting properties</b>	
No data available	

## 11.2 Information on other hazards

### Other information

No data available.

## SECTION 12: Ecological information

### 12.1 Toxicity

<b>Toxicity to fish (acute)</b>
No data available
<b>Toxicity to fish (chronic)</b>
No data available
<b>Toxicity to Daphnia (acute)</b>
No data available
<b>Toxicity to Daphnia (chronic)</b>
No data available
<b>Toxicity to algae (acute)</b>
No data available
<b>Toxicity to algae (chronic)</b>
No data available
<b>Bacteria toxicity</b>
No data available

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

### 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

Trade name: Buffer MVL

Current version : 1.0.1, issued: 06.08.2025

Replaced version: 1.0.0, issued: 02.02.2021

Region: GB

Results of PBT and vPvB assessment	
Product Name	
Buffer MVL	
PBT assessment	The product is not considered to be a PBT.
vPvB assessment	The product is not considered to be a vPvB.

**12.6 Endocrine disrupting properties**

No data available.

**12.7 Other adverse effects**

No data available.

**12.8 Other information**

Other information
Do not discharge product unmonitored into the environment.

**SECTION 13: Disposal considerations**

**13.1 Waste treatment methods**

**Product**

Disposal of the product should be carried out in accordance with all applicable regulations following consultation with the responsible local authority and the disposal company in an authorised and suitable disposal facility. Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

**Packaging**

Residues must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.

**SECTION 14: Transport information**

**14.1 UN number or ID number**

ADR/RID/ADN	UN1760
IMDG	UN1760
ICAO-TI / IATA	UN1760

**14.2 UN proper shipping name**

ADR/RID/ADN	CORROSIVE LIQUID, N.O.S.
Technical name	GUANIDINIUMTHIOCYANATE
IMDG	CORROSIVE LIQUID, N.O.S.
Technical name	GUANIDINIUMTHIOCYANATE
ICAO-TI / IATA	Corrosive liquid, n.o.s.
Technical name	GUANIDINIUMTHIOCYANATE

**14.3 Transport hazard class(es)**

ADR/RID/ADN - Class	8
Label	8
Classification code	C9
Tunnel restriction code	E
Hazard identification no.	80
IMDG - Class	8
Label	8
ICAO-TI / IATA - Class	8
Label	8

**14.4 Packing group**

ADR/RID/ADN	III
IMDG	III
ICAO-TI / IATA	III

Trade name: Buffer MVL

Current version : 1.0.1, issued: 06.08.2025

Replaced version: 1.0.0, issued: 02.02.2021

Region: GB

**14.5 Environmental hazards**

EmS F-A, S-B

**14.6 Special precautions for user**

No data available.

**14.7 Maritime transport in bulk according to IMO instruments**

Not relevant

**SECTION 15: Regulatory information****15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulations****Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)**

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances considered as substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

**REACH candidate list of substances of very high concern (SVHC) for authorisation**

According to available data and the information provided by preliminary suppliers, the product does not contain substances that are considered substances meeting the criteria for inclusion in annex XIV (List of Substances Subject to Authorisation) as laid down in Article 57 and article 59 of REACH (EC) 1907/2006.

**Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, MIXTURES AND ARTICLES**

The product is considered being subject to REACH regulation (EC) 1907/2006 annex XVII.	No 3
--	------

**Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances**

This product is not subject to Part 1 or 2 of Annex I.

**Other regulations**

Adhere to the national sanitary and occupational safety regulations when using this product.

**15.2 Chemical safety assessment**

A chemical safety assessment has not been carried out for this mixture.

**SECTION 16: Other information****Sources of key data used to compile the data sheet:**

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

Directives 2000/39/EC, 2006/15/EC, 2009/161/EU, (EU) 2017/164.

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding section.

**Full text of the H- and EUH- phrases drawn up in sections 2 and 3 (provided not already drawn up in these sections)**

H302	Harmful if swallowed.
H312	Harmful in contact with skin.
H318	Causes serious eye damage.
H332	Harmful if inhaled.

**Creation of the safety data sheet**

UMCO GmbH - D-21107 Hamburg, Georg-Wilhelm-Strasse 187, Tel.: +49(40)555 546 300, Fax: +49(40)555 546 357, e-mail: umco@umco.de

This information is based on our present knowledge and experience.

The safety data sheet describes products with a view to safety requirements.

It does not however, constitute a guarantee for any specific product properties and shall not establish a legally valid contractual relationship.

# Safety data sheet in accordance



with 1907/2006/EC

---

**Trade name:** Buffer MVL

**Current version :** 1.0.1, issued: 06.08.2025

**Replaced version:** 1.0.0, issued: 02.02.2021

**Region:** GB

---

Alterations/supplements:

Alterations to the previous edition are marked in the left-hand margin.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO GmbH.

Prod-ID 773106